CE  
1434

IUP 001-01-L-EC R06-A

SILIMED 

IMPLANTE MAMÁRIO - GEL DE SILICONE - SUPERFÍCIE LISA

PT

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1ra copia para la Paciente/ **2da copia - para doctor**

## INSTRUCTIONS ON USE / PRODUCT INFORMATION DESCRIPTION

REVIEW DATE: 29/01/2018

### DESCRIPTION

#### MAMMARY IMPLANT - SILICONE GEL - SMOOTH SURFACE

Made of a single silicone elastomer shell, mechanically resistant, with a smooth surface. It contains a defined volume of high performance transparent silicone gel (HSC) <sup>1</sup>, developed in a way that its shape and consistency provide a natural breast appearance.

The shells with low bleed treatment are a contribution to the reduction of the transudation of silicone by osmosis, by the barrier of the elastomer Shell, where part of the thickness of the Shell is formed by an inner layer of elastomer especially impermeable to the silicone oils.

It is presented in a set of types, exposed in the commercial catalog and related below, which is the result of clinical experience obtained with the medical class.

Design	REF	Basis	Profile	Projection	XXX = [VOL] (cm <sup>3</sup> )
	10512-XXX	Round	Spherical	Moderate	080; 100; 120; 140; 160; 180; 200; 220; 240; 260; 280; 300; 330; 360; 390; 420; 450; 500; 550; 600; 650; 700
	10521-XXX	Round	Spherical	High	095; 115; 135; 155; 175; 195; 215; 235; 255; 275; 295; 315; 335; 355; 385; 415; 455; 495; 545; 605; 655; 695
	10521-XXX EA	Round	Spherical	Extra High	115; 135; 155; 175; 195; 215; 235; 255; 275; 295; 315; 335; 355; 385; 415; 455; 495; 545; 605; 655; 695
	S10621 – XXX MP	Round	Spherical	Moderate	95,115,135,155,175,195,215,235,255,285,305,325,355,385,415,435,455,485,505,525,540,575,590,625,650,695
	S10621 – XXX HP	Round	Spherical	High	190,205,225,240,255,280,300,330,350,385,415,440,470,505,535,565,590,620,650,700
	S10721 – XXX MP	Round	Spherical	Moderate	95,115,135,155,175,195,215,235,255,285,305,325,355,385,415,435,455,485,505,525,540,575,590,625,650,695
	S10721 – XXX HP	Round	Spherical	High	190,205,225,240,255,280,300,330,350,385,415,440,470,505,535,565,590,620,650,700

#### <sup>1</sup> High Strength Silicone

Note 1: The difference between the models is in product dimensions and / or volumes.

It is up to the physician to select the appropriate size and type of fixation to meet the clinical and aesthetic requirements of each case.

They are all supplied sterile and intended for single use. The package should only be opened inside the surgical center.

**CONSTITUENT MATERIALS** - Raw materials are medical grade and biocompatible. The environmental conditions of manufacturing and production techniques are controlled by a strict quality management system.

- Shell: Elastomer, composition of Polydimethylsiloxane and Dimethyl Fluoro Silicone Copolymer, catalyzed by platinum compound.

- Filler: Silicone gel (HSC) <sup>1</sup>, catalyzed by a platinum complex.

The silicone may exhibit a color variation, ranging from opaque gray to translucent yellow.

**PRESENTATION** - Supplied sterile and pyrogen-free, in double packaging, packaged in sealed external carton, where the documents related to the product are found.

The label on the outer carton contains the following information, among others: product description, reference number, number of products per package, serial number, expiration date, size and / or volume of the product. The surface of the product is identified by the color used on its label, as follows:

- Smooth Surface: Green

The symbols on the label are described in the form attached to this usage instruction.

**STERILIZATION** - One of the following sterilization methods is used:

- Dry Heat or Ethylene Oxide

The sterilization process and its expiration date are indicated on the labels of the packaging, taking into account the requirements of the country to which the product is intended and which has granted the registration for marketing. Each sterilization lot receives your individual confirmation.

#### **Prohibited reprocessing.**

**INTENDED USE / PURPOSE OF USE** - Mammary Implants

- SILIMED Silicone Gel achieve the intended purpose in the following ways:

- In case of breast augmentation: increase in the volume and circumference of one or both breasts;

- In the case of breast reconstruction: return of the volume and circumference of one or both breasts removed after mastectomy or trauma, or which did not develop due to any medical abnormality.

The SILIMED mammary implant is a single use product, implantable through invasive surgery and its medical grade silicone elastomer composition and polyurethane foam coating, when applicable, come into contact with the tissues of the mammary region, aiming at a prolonged use of the implant.

**INDICATION OF USE** - Mammary Implants - SILIMED Silicone Gel - are indicated for women for the following procedures:

- Breast enlargement. Breast augmentation includes the primary increase to increase breast size, as well as revision surgery to correct or improve the outcome of a primary breast augmentation surgery.

- Breast reconstruction. Breast reconstruction includes the

primary reconstruction to replace breast tissue that was removed due to cancer or trauma or that failed to develop properly due to a serious breast abnormality. Breast reconstruction also includes revision surgery to correct or improve the results of a primary breast reconstruction surgery.

**CONTRAINDICATIONS** - The use of mammary implants or the implantation procedure is contraindicated in the following situations:

- Patients with infection in any part of the body;

- Patients with cancer or pre-malignant condition not adequately treated;

- Pregnant or breastfeeding patients;

- Patients with Atopy;

- Patients with suppressed immune system;

- Patients with local inflammation.

**PRECAUTIONS** - In the preoperative evaluation of his patients, the physician should take into account: recent history of tumors in the implant region, in particular recurrent cancer or metastasis; autoimmune diseases; pathologies and treatments that affect the immune system or blood coagulation; history of severe allergy; diabetes and cardiovascular problems.

#### **HANDLING AND USE INSTRUCTIONS**

- Opening of the package:

1 - Make sure that the outer plastic coating has not been opened;

2 - Remove the double packaging from the inside of the sealed carton. This should be carefully examined before use in the operating room. The product that has its packaging in any way violated should not be used;

3 - Separate the documents that accompany the product;

4 - Stick the adhesive labels with the product data to the records of the hospital, the doctor and also the dossier that will be delivered to the patient, as indicated by the symbols;

5 - Open the outer blister having access to the sterile inner blister containing the product. Be careful not to contaminate it on the outside of the first blister;

6 - Open the sterile inner blister in the surgical field.

The high dielectric value of silicon can generate electrostatic charges, which are responsible for attracting particles in the environment, such as plush in general and talc, to name a few. External contaminants that adhere to the surface of the Shell can generate reactions such as those caused by foreign bodies to the body, through increased fibrosis and fluid generation. For this reason, the care taken to open the carton is very important.

**SURGICAL PROCEDURE** - Appropriate surgical techniques and procedures are necessarily the responsibility of the physician, who should evaluate them according to their own training and experience, observing up-to-date and accepted techniques to minimize the occurrence of adverse reactions. The most common techniques and complications are described in this document. In addition, SILIMED has a scientific bibliography that can be requested through its representatives or directly to the head office.

It is also your physician's responsibility to make sure the staff is properly trained in order to perform the required medical

procedures.

The physician should discuss the potential risks and benefits of preoperative mammography and ultrasound for future reference. In case of suspicious changes in these exams, we suggest the investigation before or during the inclusion of the implant.

Before the surgical procedure, the patient should be informed of all procedures and all risks involved. The application of informed preoperative consents is recommended.

The physician should carefully evaluate the type of implant, access route, dissection and implant location, taking into account the patient's anatomy and the clinical and aesthetic conditions of each case.

SILIMED implants, because they are filled with highly cohesive gel silicone, can be introduced by 3 access routes: inframammary, periareolar and transaxillary. The choice of access should be in accordance with the patient's anatomy and the surgeon's experience.

The inframammary access is technically simpler, because of the easiness of access to the subglandular, subfascial and submuscular planes. It is also easier to make the placement, perform hemostasis, insert the implant and confirm its position.

The periareolar access, in half moon, requires a dissection by marking the limit of the mammary gland or sectioning it until reaching the subglandular, subfascial and submuscular planes. Depending on the thickness of the gland and the experience of the surgeon, it may be simpler or more difficult to make the placement in the chosen plane, perform hemostasis, introduce and confirm implant placement. The creation of asymmetric incisions between the breasts can lead to an asymmetry of breast enlargement. Periareolar introduction may also be associated with a higher incidence of infection and biofilm formation due to implant contact with ruptured breast ducts during implant insertion, as well as increase the risk of interference with breastfeeding.

Transaxillary access requires the dissection of a tunnel between the axillary incision and the breast space. Therefore, the distance between the axillary access and the shop boundaries independent of the plan - subglandular, subfascial and submuscular - is greater. Some surgeons may find technical difficulty in performing hemostasis and accurate dissection of the borders of the mammary sulcus by this route. Introduction of the implant requires greater skill and care in its introduction in order to avoid trauma to the axillary skin and fracture of the implant gel due to excessive force. Incorrect positioning of the implant can lead to asymmetry of sulcus, rotation and displacement of the implant.

Regardless of the route chosen, the cutaneous incision should be appropriate to the correct visualization and dissection of the placement. We suggest that the incision be between 5 and 7 cm, depending on the type of access chosen and size of the implant. The disproportion between these aspects can lead to excessive force causing burns at the edges of the incision due to the friction or fracture of the implant gel. The application of (excessive) pressure during implantation should be avoided.

Regarding the positioning, three allocations can be used for the implant: retroglandular, retrofascial and retromuscular. They all have advantages and disadvantages. These should be exposed and explained to the patient prior to implantation. As a general rule, the pre-muscle pectoral (retroglandular and retrofascial) pathways preserve the pectoral muscles and have a lower incidence of postoperative pain. However, in patients submitted to inclusion of large implants or with small breasts, the implants, especially their edges, can be palpable and visible. The retromuscular plane has thicker coverage, making it difficult to palpate and visualize the implant, but it is more painful, it impairs muscle function and can lead to poor implant placement with lateral and superior displacement.

Once the implants are inserted, the positioning of the implants must be confirmed.

It is known that the greatest immediate perioperative risks are: hematoma, seroma and infection. Thus, correct asepsis, hemostasis and implant handling are recommended in order to decrease these complications. The use of drains is possible according to the preference of the surgeon and / or need of the patient.

The basic postoperative recommendations include the use of the bra appropriate to the procedure performed, avoiding activities that may favor pain, bleeding or displacement of the implants and avoid smoking.

For postoperative follow-up, the evaluation of the patient and the implant should be made through anamnesis and physical examination, in order to detect possible changes in the implant. Imaging examinations such as ultrasound and nuclear magnetic resonance are the best so far to evaluate the integrity of the implants. We suggest performing an annual ultrasonography to evaluate the implants.

The doctor should always have a spare product available during the surgical procedures.

We suggest annual clinical follow-up.

**IMPLANT REMOVAL INSTRUCTIONS** - Removal of breast implants in general represents a simple procedure. However, possible previous complications can add some complexity to surgery.

The access route for removal of the implant should be the most appropriate at the discretion of the physician, based on previous examinations. After skin incision and dissection of the tissue planes, the plane used for implant placement should be accessed. Proceed with careful removal of the implant. Then close the planes in layers and end the surgical procedure.

It is up to the physician to choose whether or not to remove the fibrotic capsule, evaluating the impact of this withdrawal on the patient's morbidity in the postoperative period and also on possible interferences in the interpretation of mammographic exams.

In case of definitive removal, without subsequent replacement of the implant, the physician should also evaluate the esthetic result that the extraction of the capsule will cause.

In the case of removal of the implant whose Shell has suffered loss of integrity and the gel has been in contact with the tissues

of the patient, the doctor must take care to remove all the gel content and also proceed with the removal of the fibrotic capsule.

Meticulous hemostasis during surgery is recommended in order to decrease the chances of postoperative hematoma and seroma. Excessive and persistent bleeding should be monitored throughout the procedure.

The damaged product must be returned to the manufacturer for disposal.

**IMPORTANT RECOMMENDATIONS / WARNINGS** - The Mammary Implant can only be purchased by doctors or under their prescription.

It is supplied sterile and, in order to maintain its aseptic condition in the surgical act, asepsis and cleaning are required under the conditions of its use.

The implant and its packaging must be intact, otherwise the implant should not be used. It can only be used if found intact in its original form of manufacture, that is, without any change in its original characteristics.

Eventually, the implant may have blisters or a certain whitening, resulting from sterilization, which does not prevent its use.

The surface of the implant should not be contaminated with talc, dust or oils. Before contacting the implant, the medical staff should exchange the sterile surgical gloves in use with new pairs, and then wash them with sterile and apyrogenic saline solution.

Solutions containing iodine cannot come in contact with the implant.

Marker solutions, such as methylene blue, used in the medical-surgical procedure must be properly sterilized to avoid biological contamination.

Any perforating surgical instrument to be used in the region near the implant should be handled with care so that it is not reached. Any accidental bore, cut or even scratch on the implant Shell will invalidate it for surgical use. It is not possible to repair damaged product. It must be replaced by another.

The physician should instruct his patient not to massage and / or vigorously manipulate the implant or to engage in intense physical activity for a specified period of time.

The patient should be instructed to follow the most up-to-date medical recommendations regarding breast exams and to distinguish the implant from the breast tissue for the best self-examination result.

Although there is no definitive data yet, the route of scientific literature suggests that there is no direct cause and effect correlation between the placement of breast implants for aesthetic increase and frequency of suicide, and the increase in this frequency is possibly correlated with intrinsic psychiatric morbidity of these patients. Therefore, it is highly recommended that the psychological stability of the patient be carefully evaluated by the physician before and after placement of Breast Implants.

**The Breast Implant is intended for single use only. Prohibited**

**reprocessing.**

**By law, it cannot be reused or re-sterilized, as it damages the product and compromises its performance and safety.**

**The damaged product must be returned to the manufacturer for disposal.**

#### **POTENTIAL COMPLICATIONS AND POSSIBLE SOLUTIONS**

- Any patient undergoing surgery to introduce a foreign element into the body is subject to possible complications. SILIMED delegates to physicians the responsibility to clarify to patients the possibility of new surgery for removal or replacement of the implant, as well as eventual occurrence of adverse effects after its implantation and appropriate treatment. At the same time, it invites the medical profession to communicate further findings in this regard. SILIMED also delegates to the physician to adopt the best solutions available for the complications of surgery. Among these, it is noteworthy the meticulous hemostasis during the procedure, antibiotic therapy and removal or replacement of the implant.

**EXPECTED COMPLICATIONS** - For mammary implant, the potential complications of surgery, reported in the current medical-scientific literature, or directly to SILIMED, accompanied through the process of product risk management - post-production phase, are, as follows:

**THORACIC WALL SCAFFOLDING** - Few cases have been reported related to tissue expanders and, occasionally, to breast implants placed in a retromuscular position. It is associated with the surgical technique used. It can cause pain and injury to organs near the ribs and muscles of the thoracic region.

**ALERGENICITY (RASH)** - May occur in the recent postoperative period and is mainly characterized by tissue irritation, redness and rash. It is easily treatable and of transient duration.

**ALTERATION OF BREAST AND / OR NIPPLE SENSITIVITY** - Uncommon. Corresponds to the increase (hyperesthesia) or decrease (hypoesthesia) of sensitivity of the implanted region due to the surgical incision. It is usually a reversible process over time, but it can last for two years or more.

**TISSUE ATROPHY** - Associated with large volume breast implants, when inserted in stores disproportionate to the volume of these implants. Mechanical forces press the parenchyma and subcutaneous tissue, which can cause atrophy of these tissues over time.

**CALCIFICATION** - Characterized by the presence of deposits of insoluble salts of calcium (and / or magnesium), especially calcium carbonate and phosphate in the fibrotic capsule. It can cause pain and hardening of the healing tissue. Of unknown origin, it was observed in the tissues neighboring the implantation ones, sometimes, requiring the removal of the implant.

**CYSTS / NODULES** - Nodules and breast cysts are changes that occur in the breast tissue. Breast lumps are a bulge of the breast. They can usually be caused by fibrocystic changes, fibroadenomas, cysts and cancer. They can also occur due to implant rupture if extravasation occurs and the silicone gel migrates to lymph nodes or other parts of the body. The cyst is characterized as a closed sac / cavity filled with fluid coated by an epithelium.

**IATROGENIC COMPLICATIONS** - Certain situations such as folds and perforations in the breast implant, among others, can be caused by the use of inadequate surgical technique.

**CAPSULE CONTRACTURE / FIBROSIS** - Occurs when the healing tissue or fibrotic capsule becomes thicker and retractable and compresses the implant. The fibrotic capsule forms around the implant due to the reaction of the organism to the presence of a foreign body. Factors such as infection, hematoma and seroma may favor the appearance of capsular contracture. The likelihood of occurrence decreases with meticulous hemostasis of the incision and careful handling of the implant. SILIMED advises against an external (closed) capsulotomy procedure, which may lead to rupture of the implant.

**SUTURE DEHISCENCE** - Suture dehiscence is the separation of the layers of the surgical wound. It may be partial and superficial or complete when there is complete disruption of all layers of the tissue. Normally, this effect starts with a small seroma. In most cases it is a problem related to the surgical technique and / or use of a larger implant, more than the adequate.

**IMPLANT DISPLACEMENT / BAD POSITIONING** - It can happen with patient discomfort and distortion of the desirable shape, a problem related to the operative technique, ineffective formation of the fibrotic capsule or trauma occurrence. The physician should evaluate the need for removal of the implant and / or its probable replacement.

**PAIN** - Some patients may experience pain in the postoperative period, which should be investigated immediately. It can vary from mild to severe and from short to long lasting. It can result from infection, inflammation, seroma, calcification and / or use of inadequate surgical technique.

**WRINKLE** - Ripple of the surface of the breast implant can occur, which can be seen and / or felt.

**EROSION** - After implantation, erosion of tissue near the implant can occur, caused by folds in the device and other deformities that lead to exaggerated tissue compression.

**EXTRUSION** - Unstable or compromised tissue coverage, excessive stress or trauma at the implant site, and / or disruption of the healing process can lead to implant extrusion. It has been reported as an uncommon complication in both patients who have undergone augmentation and reconstruction surgery.

**GALACTORRHEA** - May be the result of hormonal disorders or the action of medications. It causes discomfort and embarrassment, affecting the emotional state of the patient.

**GRANULOMA** - It is associated with the formation of granular tissue due to the presence of various cell types such as macrophages, lymphocytes and fibroblasts. It can lead to a chronic inflammatory lesion, tumor and consequent surgical intervention to remove the implant.

**HEMATOMA** - Extravasation of blood and its confinement in an organ or tissue (hematoma) can occur in any type of surgery. In its occurrence, it must be carefully evaluated to define the conduct to be adopted. At risk of being a precursor of infection and increased fibrosis, if not resorbed, the hematoma should be removed. The main preventive measure is perfect hemostasis during surgery and postoperative rest.

**INFECTION (BIOLOGICAL CONTAMINATION)** - Uncommon and difficult due to implantation, if the recommendations of use and asepsis associated with implantation surgery are observed. It includes the possibility of occurrence of secondary necrosis and extrusion of the implant. Local culture is recommended and, in sequence, the use of antibiotics, drainage and, eventually, removal of the implant.

**INTERFERENCE WITH MAMMOGRAPHY; ULTRASOUND AND MAGNETIC RESONANCE IMAGING** - Interference occurs when the presence of implants (or artifacts) causes an erroneous reading of the result, for example: false positive diagnosis of rupture, resulting in unnecessary removals. Breast implants can interfere with the diagnosis of mammography breast cancer by covering up some suspicious spot on the breast tissue, masking the results and delaying the onset of anti-cancer treatment.

**BREASTFEEDING INTERFERENCE** - The presence of silicone in breast milk was investigated in several studies, and no significant levels of the substance were found in mothers with silicone breast implants when compared to women without implants. There is, however, a possible interference with breastfeeding (reduction in milk production) in women who underwent mastoplasty with breast implants. The frequency of breastfeeding interference is greater in patients in whom the periareolar incision of the breast implant was used.

**NECROSIS** - Result of exaggerated skin tension caused by the presence of the implant or by surgical trauma. The physician should evaluate the need for removal of the implant and / or its probable replacement.

**LOSS OF MECHANICAL INTEGRITY / PRODUCT FAILURE** - In case of loss of mechanical integrity of the breast implant, the patient may undergo surgery for replacement, depending on the degree of the hazard. Among these situations, it is possible to mention the implant detachment, the spontaneous rupture of its Shell, among others.

**PNEUMOTHORAX** - This is the collapse of the lung due to a sudden change of pressure inside the thoracic cavity. In the case of breast implants, this can occur when the implantation shop is formed due to surgical trauma, mainly in the retromuscular space. It has the following consequences: tachycardia, chest pain, cough, shortness of breath, cyanosis, agitation, among others.

**INFLAMMATORY RESPONSE** - The presence of a foreign element in the body triggers an inflammatory process, which may be more or less intense. Severe inflammation may also be due to trauma or inadequate procedures. Causes pain, flushing, local heat and swelling. This effect can be treated pharmacologically.

**UNSATISFACTORY AESTHETIC RESULT** - It is associated with patient and / or physician dissatisfaction with the result of breast implant placement surgery. Some possible unsatisfactory aesthetic results are: ptosis, sagging, asymmetry, folds, bends, implant wrinkling, stretch marks, adhesion and hypertrophic scars. The physician should evaluate the best measure to be taken, and there is no need for surgical intervention for implant replacement.

**RUPTURE** - It can be generated by the natural wear of the implant by implantation time, by puncture-cutting instrument damage, by excessive effort during surgery, by traumas or extreme physical pressure, by compression during mammography or by severe capsular contracture. Most occur during the surgical procedure. In case of suspected rupture of the implant in an asymptomatic patient ("silent rupture"), a careful evaluation should be made to confirm the situation and decide the procedures to be taken, including the possibility of implant replacement.

**SEROMA** - Accompanied by edema and pain, and it can be

relieved with rest, immobilization and application of compresses, eventually drainage and even removal of the implant. The seroma is not always sharp. There may be a minimal, continuous, undetected seroma that appears until months later as a "late seroma." There are also "late seromas" of unknown origin.

**SYMPTOMS ASSOCIATED WITH MACROMASTÍA** - They are associated with breasts of large volumes. Some of the symptoms are pain in the following areas: back, breasts, hands (including numbness), neck and head; skin irritation and marks on the shoulders. The volume of the breast implant should be considered respecting the amount of breast tissue present, the shape of the breasts and the patient's desire. The methodologies used to estimate if the volume of breast implants is adequate to patients' needs and expectations are of great importance in the choice of breast implant.

**GEL BLEEDING** - It is associated with a silicone gel filled breast implant. It is the process of diffusion and migration of the gel through the entire Shell (transudation) or by rupture of the implant into adjacent tissues. It may cause local complications such as: reduction in size, asymmetry, alteration of breast shape, pain, capsular contracture or formation of nodules and granulomas. **SILIMED** high-performance silicone gel minimizes the occurrence of this adverse effect.

**TROMBOSIS** - Uncommon, usually associated with the surgical procedure or patients with obesity, neoplasm, varicose veins, traumatic or iatrogenic venous lesions, bedridden patients or immobilized by a systemic or local disease in the limb, patients with heart disease, extrinsic compression of the vein, pregnancy and puerperium, use of hormonal contraceptives and hematological diseases. It may lead to local complications and removal of the implant may be necessary.

Other complications: Emotional alteration, Discomfort, Numbness, Edema, Erythema, Hardening of the breast, Hemorrhage, Swollen lymph nodes (axillary lymph nodes) also called Lymphadenopathy, Tissue Infiltration, Irritation, Tissue Lesion, Palpability, Pruritus, Flushing, Skin Tightness, Implant visibility, Asymmetry, Ptosis, Late cicatrization, Upper pole filling.

Some complications have been reported in the medical-scientific literature, but to date there is no scientific evidence of a possible cause and effect relationship between these events and silicone implants. Some examples are:

**CANCER** - According to the current medical-scientific literature, women undergoing silicone implant surgery are at no greater risk of developing cancer.

**ANAPLASTIC LARGE CELL LYMPHOMA (ALCL)** - Not a breast cancer, it is a rare type of Non-Hodgkin's Lymphoma. Women with breast implants may have a small risk of developing ALCL in the fluid (seroma) or in the fibrous capsule adjacent to the implant. Based on the information available, it is not possible to confirm with statistical certainty that breast implants cause ALCL.

**SILIMED** supports FDA recommendations to physicians and patients advising physicians to remain diligent in monitoring any symptoms of ALCL in patients with breast implants regardless of their rare occurrence by encouraging the health care provider to report all confirmed cases of ALCL in women with breast implants to the Regulatory Agency of their country.

In addition, prophylactic removal of breast implants in patients

with no symptoms or other abnormalities is not recommended. The healthcare provider should consider the possibility of ALCL if the patient has late healing, persistent fluid around the implant, and if there are changes in the region around the implant, such as pain and swelling. In cases of seroma, it is recommended to send the freshly collected fluid for pathological analysis in order to rule out the possibility of ALCL. Samples of tissue around the implant should also be analyzed. If diagnosed, an individual treatment will be performed, because, due to the occurrence of a small number of cases in the world and a variety of treatments, there is no single defined treatment. Women should monitor their breast implants and contact their physicians if they notice any changes, and those who are considering breast implant surgery should discuss the risks and benefits with their physicians.

**IMMUNOLOGICAL RESPONSE** - Studies to date have found no evidence that the placement of silicone implants causes connective or autoimmune diseases. Current literature suggests a possible link between implant rupture and fibromyalgia and connective tissue diseases. Further studies are needed.

**TERATOGENICITY OR OTHER REPRODUCTIVE EFFECTS** - It is the ability to cause congenital malformations in the fetus. The review of the medical-scientific literature indicates that studies have shown no evidence of teratogenicity or other reproductive effects associated with breast implants.

**EFFECTS ON CHILDREN WHO ARE BORN FROM MOTHERS WITH MAMMARY IMPLANTS** - Recent studies have shown that the amount of silicone in breast milk found in women who have silicone implants is not significantly different from women who do not have silicone implants. However, the review of the medical-scientific literature does not demonstrate evidence that silicone gel breast implants can cause any harmful effects to the offspring of women with implants.

**SUICIDE** - Epidemiological studies indicate that women who underwent augmentation mammoplasty have a suicide rate greater than the general population. Although there is no definitive data yet, the scientific literature suggests that there is no direct cause and effect correlation between the placement of breast implants for aesthetic increase and suicide frequency, and the increase in this frequency is possibly correlated with the intrinsic psychiatric morbidity of these patients. Therefore, it is highly recommended that the psychological stability of the patient be carefully evaluated by the physician before and after placement of Breast Implants.

**DURABILITY** - Like all breast implants, **SILIMED** Mammary Implants have limited useful life. This lifetime has not been determined by the scientific community yet. **SILIMED** establishes a mean period of 10 years as the expected useful life, but such parameter may be altered whenever some reason to justify it is assessed. The implant may need to be removed or replaced, which may imply in a revision surgery. An annual clinical follow-up is recommended for these patients.

**INFORMATION FOR SILIMED PATIENT** - The following information is directed to all patients who wish to undergo breast augmentation or reconstructive surgery using silicone breast implants. Information on the main questions about the subject, including potential complications, benefits and risks, activities that may damage the implant, the possible need for removal / replacement of devices and the expected durability of the breast implant.

The physician responsible for the procedure should carefully read this section and be responsible for all patient questions prior to the SILIMED mammary implant procedure.

Annual follow-up of these patients is recommended. In case of any suspected complications listed in this information or any abnormal symptoms, contact your doctor.

**EXPECTED BENEFITS** - A mammary implant surgery can bring great benefit by providing successful breast reconstruction, augmentation and surgical revision. Doing breast reconstruction / augmentation also results in improvements in quality of life. In addition, studies have shown that breast reconstruction with breast implants has been an aid in breast cancer recovery, as well as reducing emotional stress by helping the body to return to a more natural appearance, as opposed to not having surgery or using external prosthesis.

**PATIENT CLARIFICATION AND CONSENT** - Considering the risks inherent in a surgery, with or without the use of implants, and the possible complications related to them, SILIMED relies on physicians to clarify their patients about the existing risk-benefit balance.

Each patient should receive SILIMED Patient Information during their initial consultation to allow sufficient time prior to surgery to properly read and understand the important risk information, follow-up recommendations, and benefits associated with breast implant placement surgeries. To document the patient's clarification successfully, the Term of Clarification for the use of the SILIMED Mammary Implant must be signed by the patient and the physician and attached to the patient's medical record.

**INTERFERENCE WITH MAMMOGRAPHY; ULTRASOUND AND MAGNETIC RESONANCE IMAGING** -Interference occurs when the presence of implants (or artifacts) causes a false reading of the result, for example: false positive diagnosis of rupture, resulting in unnecessary withdrawals.

SILIMED points out that the implant may interfere with the quality of mammography. Therefore, the patient should be advised and instructed to request the services of professionals who have experience with mammography techniques with implant, informing the radiologist the need to adapt the mammographic compression and not forgetting to show the "SILIMED PRODUCT TRACEABILITY CARD - FOR PATIENT".

Other methods, such as ultrasound and magnetic resonance imaging may be useful in conjunction with mammography because they do not require compression and allow examination from any angle.

According to recent studies, silicone breast implants reduce the image quality of echocardiographic visualizations. More data is needed, however. SILIMED advises its patients to inform the technician responsible for the examination of the existence of silicone implants.

**CAUTIONS** - Within the first month after surgery, some activities can impair the implant and should be avoided, such as:

- Exposure to the sun;
- Sudden movements;
- Sports in general.

You should be aware that normal stress or trauma applied at the surgery placement can cause implant extrusion.

The use of proper bra and not doing violent exercises are minimal recommendations to be followed in the postoperative period.

The patients should inform the medical or pharmaceutical specialist that they are carrying silicone implants before using topical medicine in the breast region.

The patients must continue to consult with the medical specialist to continue routine monitoring for breast cancer.

The patient should inform the physician about the presence of the implant if any breast surgery is scheduled.

If you suspect of any of the complications listed in this information or any abnormal symptoms, especially in case of trauma or compression (e.g. excessive massage of the breast through some sports or wearing a seat belt), contact your doctor.

#### **IMPORTANT FACTORS TO CONSIDER** -

- 1.1. If you are going to undergo breast augmentation or reconstruction, be aware that breast implants are not considered lifelong devices and that breast implantation may not be a single surgery. You may need one or more additional surgeries and have to schedule doctor visits throughout your life. You may also need surgery to remove the implant, with or without replacement.
- 2.2. Many of the changes in your breast after implantation cannot be undone. If you subsequently choose to have the implant (s) removed, you may experience unacceptable chest wall deformity, unacceptable crimping, unacceptable wrinkling, loss of breast tissue or other changes in the breast contour.
- 3.3. With breast implants, routine mammography screening will be more difficult and you will need additional imaging, which means more radiation and more time.
- 4.4. Every SILIMED Mammary Implant comes with a "Silimed Product Traceability Card - For Patient" which contains all the information regarding the surgery. This card is for your safety, you must carry this card to facilitate medical care in an emergency.

**STORAGE/CONSERVATION/TRANSPORTATION** - SILIMED products and packaging are resistant and, if stored at ambient temperature and in accordance with the guidelines printed on the packaging, are not damaged. Because they are sterile medical products, it is recommended that they not be transported or stored together with other types of material that could cause physical or chemical damages to the product, which would prevent its use.

**WARRANTY** - SILIMED will replace any of its line products that presents a provable manufacturing defect, provided it is returned duly identified by the person who purchased it and after analysis at the SILIMED laboratory.

The adverse effects presented here, the inadequate surgical use and color variation mentioned in the item "Constituent Materials", are not considered defects of the product.

The guarantee of SILIMED products does not cover the simple decision of the patient or the doctor of the product exchange.

**TRACEABILITY** - SILIMED keeps accurate records, unit by unit, of raw materials, manufacturing steps, environmental conditions, operating conditions and quality tests, assigning a unique number to each product, which fully identifies it at any time. In this sense, it is absolutely necessary that any claim be forwarded with the SERIAL NUMBER (SN) or BATCH NUMBER, when applicable, of the claimed unit, which appears on the packaging.

SILIMED provides adhesive labels with the product data that must be attached to the hospital and doctor records as indicated by the symbols.

Also provided are: Patient Registration Form, to be signed by the patient and returned to SILIMED by the physician, and Traceability Card of SILIMED products, to be delivered by the doctor to the patient. These documents contain product data, which enables traceability.

**PRODUCT RETURN** - Products contaminated with biological material must be properly decontaminated by the physician or responsible health unit prior to shipment, and returned to SILIMED following applicable legislation.

No product must be returned without prior authorization from SILIMED.

Dear Doctor: In order to facilitate understanding when handling our products, in this printed notice you will find international symbols, together with those developed at SILIMED, that appear on our labels, packaging, instructions for use and/or catalogues, as well as their meanings.

	Single use		Do not use if package is damaged
	Lot number		Volume of the product
	Date of manufacture		Consult instructions for use
	Ethyleneoxide sterilization		Keep away from sunlight
	Steam or dry-heat sterilization		Dimensions of product
	Manufacturer		CE mark
	Use by		Quantity of product
	Serial number		ANVISA registration number
	Sterile		Label for tracing to the doctor's dossier
	Reference number		Verified by the QC

	Caution, consult accompanying documents		INMETRO seal for identification of conformity
	Authorised representative in the european community		Label for tracing to the hospital dossier
	Label for tracing the patient's dossier		Medical device
	Keep dry		Temperature limit
	This way up		Humidity limit
	Right Side		Unique Device Identifier
	Left Side		Portuguese
	Spanish		Italian
	Russian		English

## INFORMATION FOR PATIENT - SILIMED MAMMARY IMPLANTS

Silimed has developed this booklet to help you understand the benefits and risks related to breast augmentation or reconstructive surgery using silicone breast implants. After you have read this leaflet carefully, talk to your doctor before making any decisions. If you decide on the procedure, you and your doctor should sign the Clarification Statement before surgery confirming that you have read and understood the information provided about the benefits and risks of silicone breast implant surgery. The Term of Clarification for the use of the SILIMED Breast Implant is found at the end of this leaflet.

**DESCRIPTION** - SILIMED Mammmary Implants consist of a single outer elastomer membrane, filled with medical grade silicone gel; may have a smooth surface, textured or even coated with polyurethane foam. They are all supplied sterile and intended for single use.

These implants are available in various sizes and shapes to help each woman achieve the best result for her body. The doctor should guide you in choosing the correct implant to achieve the desired esthetic effect.

**EXPECTED BENEFITS** – Breast augmentation / reconstruction surgery can bring great benefit by providing reconstruction, augmentation and surgical revision of the breast with a high index of satisfaction. Patients choose primary breast augmentation surgery to increase the size and proportion of their breast(s). In addition, patients choose revision surgery (replacement of an existing breast implant) to correct or improve the results of a primary augmentation surgery. According to the medical-scientific literature, most women who underwent breast implant surgery reported high levels of satisfaction with body image, shape and size of the breasts. Additionally, after surgery, the patients reported higher self-esteem, and they felt more feminine and attractive. Breast reconstruction surgery also results in improved quality of life. In addition, studies show that breast reconstruction with breast implant has been an aid in the recovery of breast cancer.

**CLARIFICATION AND CONSENT OF THE PATIENT** - Considering the risks inherent to surgery, with or without the use of implants, and the possible complications related to them, SILIMED counts on physicians to clarify their patients about the existing risk-benefit balance.

Each patient should receive Patient Information - SILIMED Mammmary Implants during their initial consultation to allow sufficient time prior to surgery to properly read and understand important risk information, follow-up recommendations, and benefits associated with breast implant placement surgeries. To document the patient's clarification successfully, the Term of Clarification for the use of the SILIMED Breast Implant must be signed by the patient and the physician and attached to the patient's medical record.

**FORESEEN COMPLICATIONS** - Breast implant surgery offers great benefits. However, any type of surgery involves risks. There are some local complications (problems near the surgery placement) that may occur after surgery. SILIMED delegates to physicians the responsibility to inform patients about the possible occurrence of adverse effects, consequences, adequate treatment. The following are the potential risks and complications reported in the literature associated with breast implant surgery in different models and manufacturers, as well as the possible effects related to these risks.

Alteration of breast and / or nipple sensitivity, Thoracic wall scaffolding, Allergenicity / Rash, Emotional Change / Depression, Tissue Atrophy, Calcification, Cysts / nodules / tumor, Iatrogenic Complications, Capsular Contracture / Fibrosis, Suture Dehiscence, Discomfort, Implant Displacement, Pain, Numbness, Edema (Swelling), Erythema ( Hyperemia), Breast Hardening, Wrinkling, Erosion, Extrusion, Galactorrhea, Granuloma, Hematoma, Hemorrhage, Implant Displacement, Pain, Numbness, Edema (Swelling), Erythema (Hyperemia), Breast Enhancement, Wrinkling, Erosion, Extrusion, Galactorrhea, Granuloma, Hematoma, Hemorrhage, Implant Removal, Infection (Biological Contamination) and its symptoms, Lymphatic Swelling (Axillary Lymph Nodes), Tissue Infiltration, Irritation, Tissue Injury, Breastfeeding Interference, Malposition, Silicone Migration (Lymphadenopathy and Silicoma), Necrosis, Loss of Mechanical Integrity / Product Failure, Palpability, Pneumothorax, Pruritus, Local Inflammatory Response, Revision Surgery, Unsatisfactory Aesthetic Result (such as such as ptosis, sagging, asymmetry, folds, wrinkles, adhesions, stretch marks, and hypertrophic scars, Blushing, Breaking (including cases of "silent" rupture), Secretion, Seroma / fluid, Symmastia, Macromastia Associated Symptoms, Skin Tension, Gel Transudation (gel bleeding), Thrombosis, Vesicles, Implant visibility, and Top pole filling.

Other complications have been reported in the medical-scientific literature, but to date there is no scientific evidence of a possible cause and effect relationship between these events and silicone implants. Some examples are:

**CANCER** - According to the current medical-scientific literature, women undergoing silicone implant surgery are at no greater risk of developing cancer.

**ANAPLASTIC LARGE CELL LYMPHOMA (ALCL)** - Not a breast cancer, it is a rare type of Non-Hodgkin's Lymphoma. Women with breast implants may have a small risk of developing ALCL in the fluid (seroma) or in the fibrous capsule adjacent to the implant. Based on the information available, it is not possible to confirm with statistical certainty that breast implants cause ALCL.

SILIMED supports FDA recommendations to physicians and patients advising physicians to remain diligent in monitoring any symptoms of ALCL in patients with breast implants regardless of their rare occurrence by encouraging the health care provider to report all confirmed cases of ALCL in women with breast implants to the Regulatory Agency of their country.

In addition, prophylactic removal of breast implants in patients with no symptoms or other abnormalities is not recommended. The physician should consider the possibility of ALCL if the patient has late healing, persistent fluid around the implant, and if there are changes in the region around the implant, such as pain and swelling.

In cases of seroma, it is recommended to send the freshly collected fluid for pathological analysis in order to rule out the possibility of ALCL. Samples of tissue around the implant should also be analyzed. If diagnosed, an individual treatment will be performed, because due to the occurrence of a small number of cases in the world and a variety of treatments, there is no single defined treatment. Women should monitor their breast implants and contact their physicians if they notice any changes, and those who are considering breast implant surgery should discuss the risks and benefits with their physicians.

**IMMUNOLOGICAL RESPONSE** - Studies to date have found no evidence that the placement of silicone implants causes connective or autoimmune diseases. Current literature suggests a possible link between implant rupture and fibromyalgia and connective tissue diseases. Further studies are needed.

**TERATOGENICITY OR OTHER REPRODUCTIVE EFFECTS**- It is the ability to cause congenital malformations in the fetus. The review of the medical-scientific literature indicates that studies have shown no evidence of teratogenicity or other reproductive effects associated with breast implants.

**EFFECTS ON CHILDREN WHO ARE BORN FROM MOTHERS WITH MAMMARY IMPLANTS**- Recent studies have shown that the amount of silicone in breast milk found in women who have silicone implants is not significantly different from women who do not have silicone implants. However, the review of the medical-scientific literature does not demonstrate evidence that silicone gel breast implants can cause any harmful effects to the offspring of women with implants.

**SUICIDE** - Epidemiological studies indicate that women who underwent augmentation mammoplasty have a suicide rate greater than the general population. Although there is no definitive data yet, the scientific literature suggests that there is no direct cause and effect correlation between the placement of breast implants for aesthetic increase and suicide frequency, and the increase in this frequency is possibly correlated with the intrinsic psychiatric morbidity of these patients. Therefore, it is highly recommended that the psychological stability of the patient be carefully evaluated by the physician before and after placement of Breast Implants.

**INTERFERENCE IN BREASTFEEDING** – The presence of silicone in the mother’s milk has been researched in several studies without significant levels of the substance being found in mothers who bear silicone mammary implants when compared with women without implants. There is, however, a possible interference in breast-feeding (reduced milk production) in women who have undergone surgery with mammary implants. The frequency of interference in breast-feeding is higher for patients whose incision route was periareolar.

**INTERFERENCE WITH MAMMOGRAPHY; ULTRASOUND AND MAGNETIC RESONANCE IMAGING** - Interference occurs when the presence of implants (or artifacts) causes an erroneous reading of the result.

SILIMED points out that the implant may interfere with the quality of mammography. Therefore, you should be guided and instructed to request the services of professionals who have experience with mammography techniques with implants, informing the radiologist of the need to adapt the mammographic compression and not forgetting to show the “PRODUCT TRACEABILITY CARD SILIMED - FOR PATIENT “.

Other methods, such as ultrasound and magnetic resonance imaging, may be useful in conjunction with mammography because they do not require compression and allow examination at any angle.

According to recent studies, silicone breast implants reduce the image quality of echocardiographic visualizations. More data is needed, however. SILIMED advises its patients to inform the technician responsible for the examination of the existence of silicone implants.

**INTERFERENCE WITH SELF-EXAMINATION** - Breast self-examination is a prevention technique used to identify early stages of breast cancer.

The doctor should advise you on how to distinguish the implant from the breast tissue during self-examination.

You should perform the self-examination periodically for tracking lumps, swelling, hardening, or change in implant shape, which may be signs of rupture. If any of these symptoms or persistent pain occur, tell your doctor.

It is important to emphasize that breast examination done by the woman herself does not replace the physical examination performed by a qualified health professional (doctor or nurse) for this activity.

**PRECAUTIONS** - Tell your doctor if you have any of the following conditions:

- Autoimmune disease;
- Recent history of cancer;
- Weakened immune system (users of immunosuppressive drugs);
- Conditions that interfere with blood clotting;
- History of severe allergy;
- Cardiovascular diseases;
- Diabetes;

Before surgery, you should have a detailed discussion with your doctor about your medical history.

**CAUTIONS** - In the first month after surgery, some activities can impair the implant and should be avoided, such as:

- Exposure to the sun;
- Sudden movements;
- Sports in general.

You should be aware that normal stress or trauma applied at the surgery placement can cause implant extrusion.

The use of proper bra and not doing violent exercises are minimal recommendations to be followed in the postoperative period.

Ask your doctor about the activities he does not recommend in the postoperative period.

You should tell the medical or pharmaceutical specialist that you are carrying silicone implants before using topical medicine (for example, steroids) in the breast region.

You need to continue to consult with the medical specialist to keep with the monitoring routine to detect breast cancer.

You should tell the doctor about the presence of the implant if any breast surgery is scheduled.

In case of any suspected complications listed in this information or any abnormal symptoms, contact your doctor.

**DURABILITY** - Like all breast implants, SILIMED Mammary Implants have limited useful life. This lifetime has not been determined by the scientific community yet. SILIMED establishes a mean period of 10 years as the expected useful life, but such parameter may be altered whenever some reason to justify it is assessed.

The implant may need to be removed or replaced, which may imply in a revision surgery. An annual clinical follow-up is recommended for these patients.

#### **IMPORTANT FACTORS TO CONSIDER:**

1. If you are going to undergo breast augmentation or reconstruction, be aware that breast implants are not considered lifelong devices and that breast implantation may not be a single surgery. You may need one or more additional surgeries as well as medical follow-up throughout your life. You may also need surgery to remove the implant, with or without replacement.

2. Many of the changes in your breast after implantation cannot be undone. If you subsequently choose to have the implant (s) removed, you may experience unacceptable chest wall deformity, unacceptable crimping, unacceptable wrinkling, loss of breast tissue or other changes in the breast contour.

3. With breast implants, routine mammography screening will be more difficult and you will need additional imaging, which means more exposure to radiation.

4. Every SILIMED Mammary Implant comes with a "SILIMED PRODUCT TRACEABILITY CARD - FOR PATIENT", which contains all the information regarding the surgery. This card is for your safety, you must carry this card to facilitate medical care in an emergency.

This document is available on our website: [www.silimed.com.br](http://www.silimed.com.br)

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# CLARIFICATION TERM FOR THE USE OF SILIMED MAMMARY IMPLANTS

## Dear Patient,

In order to support the decision on the use of the SILIMED Mammary Implant, we will make available the Patient Information - SILIMED Mammary Implants on the company website, so that you can access important issues and clarify any doubts with your doctor before undergoing surgery.

Before and after surgery, you should consult the doctor and undergo periodic examinations to monitor their health conditions.

It is important that you be aware of the need to follow the guidelines provided by your doctor correctly to be successful in the surgery of placement of SILIMED Mammary Implant. Breast implants have a limited lifetime and may need to be removed or replaced, which may require revision surgery.

Both the physician and the anesthesiologist must have your consent to take the necessary steps to restore their clinical condition, if necessary.

The success of surgery with the MammaryImplant SILIMED depends on the commitment to follow the medical recommendations, which will allow the increase of the projection or reconstruction of the breast and without major complications. Breast surgery can improve your self-esteem and quality of life.

## Confirmation

I, \_\_\_\_\_, identity document no \_\_\_\_\_, read the information for the patient - SILIMED MAMMARY IMPLANTS made available by SILIMED, I understood the product information and surgery procedures and had the opportunity to clarify with my physician any doubts about the nature of the product, the its benefits, risks and potential compli-cations, procedures and objectives of surgery. I am aware of the possibility of adverse effects, such as those provided in the Patient Information SILIMED, and that the occurrence of adverse effects does not characterize the product defect for warranty purposes.

Patient's Signature \_\_\_\_\_ date \_\_\_\_/\_\_\_\_/\_\_\_\_ IID Number:: \_\_\_\_\_

Witness \_\_\_\_\_ date \_\_\_\_/\_\_\_\_/\_\_\_\_ ID Number: \_\_\_\_\_

Physician \_\_\_\_\_ date \_\_\_\_/\_\_\_\_/\_\_\_\_ Regional Board of Medicine: \_\_\_\_\_

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1st copy for the Patient / 2nd copy for the Physician

## Dear Patient,

In order to support the decision on the use of the SILIMED Mammary Implant, we will make available the Patient Information - SILIMED Mammary Implants on the company website, so that you can access important issues and clarify any doubts with your doctor before undergoing surgery.

Before and after surgery, you should consult the doctor and undergo periodic examinations to monitor their health conditions.

It is important that you be aware of the need to follow the guidelines provided by your doctor correctly to be successful in the surgery of placement of SILIMED Mammary Implant. Breast implants have a limited lifetime and may need to be removed or replaced, which may require revision surgery.

Both the physician and the anesthesiologist must have your consent to take the necessary steps to restore their clinical condition, if necessary.

The success of surgery with the MammaryImplant SILIMED depends on the commitment to follow the medical recommendations, which will allow the increase of the projection or reconstruction of the breast and without major complications. Breast surgery can improve your self-esteem and quality of life.

## Confirmation

I, \_\_\_\_\_, identity document no \_\_\_\_\_, read the information for the patient - SILIMED MAMMARY IMPLANTS made available by SILIMED, I understood the product information and surgery procedures and had the opportunity to clarify with my physician any doubts about the nature of the product, the its benefits, risks and potential compli-cations, procedures and objectives of surgery. I am aware of the possibility of adverse effects, such as those provided in the Patient Information SILIMED, and that the occurrence of adverse effects does not characterize the product defect for warranty purposes.

Patient's Signature \_\_\_\_\_ date \_\_\_\_/\_\_\_\_/\_\_\_\_ IID Number:: \_\_\_\_\_

Witness \_\_\_\_\_ date \_\_\_\_/\_\_\_\_/\_\_\_\_ ID Number: \_\_\_\_\_

Physician \_\_\_\_\_ date \_\_\_\_/\_\_\_\_/\_\_\_\_ Regional Board of Medicine: \_\_\_\_\_

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1st copy for the Patient / **2nd copy for the Physician**

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