IMPLANTE MAMÁRIO - GEL DE SILICONE - SUPERFÍCIE TEXTURIZADA

PT

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1ra copia para la Paciente/ 2da copia - para doctor

INSTRUCTIONS FOR USE / INFORMATION ON THE PRODUCT

Review Date: 13/07/2020

DESCRIPTION

MAMMARY IMPLANT – SILICONE GEL – TEXTURED SURFACE

It consists of a single, mechanically resistant silicone elastomer shell with a textured surface. It contains a defined volume of transparent, high performance silicone gel (HSC) ¹, developed in such a way that its shape and consistency provide a natural appearance to the breast.

The shells with low bleed treatment are a contribution to the reduction of silicone transudation by osmosis, through the barrier of the elastomeric shell, where part of the shell thickness is formed by an inner layer of elastomer especially impermeable to silicone oils.

The textured surface is effective in reducing the occurrence of capsular contracture, according to scientific evidence, allying itself with the low bleed shell in the greater ability to prevent this clinical condition.

It is presented in a set of types, which is the result of the clinical experience obtained from the medical class.

Note 1: For more information regarding the implants, see the table available on the cover of this instruction for Use.

Note 2: The Natural, Nuance and Enhance models can have a position indicator on the anterior and / or posterior sides, to facilitate the positioning of the product during the surgical procedure.

Note 3: The difference between the models is in the dimensions and / or volumes of the product.

Accessories: The Breast Implant - Silicone Gel - Textured Surface is supplied with an auxiliary glove for implantation. This accessory is not sold separately.

It is up to the physician to properly select the size and type of fixation to meet the clinical and aesthetic requirements of each case.

They are all supplied sterile and intended for single use. The packaging should only be opened inside the operating room.

CONSTITUENT MATERIALS -The raw materials are medical grade and biocompatible. The environmental conditions of manufacture and production techniques are controlled by a strict quality management system.

- Shell: Elastomer, composition of Polydimethylsiloxane and Dimethyl Fluoro Silicone Copolymer, catalyzed by a platinum compound.

- Textured surface: Formed by deposition of layers of Polydimethylsiloxane, which creates a rough surface, forming an open cell structure.

- Filling product: Silicone gel (HSC) ¹, catalyzed by a platinum complex.

Silicone can vary in color, from opaque gray to translucent yellow.

PRESENTATION - Supplied sterile and apyrogenic, in double packaging, packed in a sealed external box, where the

documents related to the product can be found.

The outer packaging label contains, among others, the following information: product description, reference number, quantity of products per package, serial number, expiration date, size and / or volume of the product. The product surface is identified by the color used on its label, as follows:

- Textured Surface: Blue

The symbols shown on the label are described at the end of this instruction for use.

STERILIZATION - One of the following sterilization processes is used:

Dry-heat or Ethylene Oxide

The sterilization process and its expiry date are indicated on the package labels. They meet the requirements of the country to which the product is destined and which granted the register for commercialization. Each sterilization lot receives individual confirmation.

Prohibited reprocessing.

INTENDED USE / PURPOSE OF USE - SILIMED Silicone Gel Mammary Implants achieve their intended purpose in the following ways:

- in case of breast enlargement: increased volume and circumference of one or both breasts;

- in the case of breast reconstruction: return of the volume and circumference of one or both breasts removed after mastectomy or mesh or not developed due to any medical abnormality.

The SILIMED Mammary Implant is a single-use product implantable through invasive surgery and its medical grade silicone elastomer composition and polyurethane foam lining, where applicable, come into contact with the tissues of the mammary region for prolonged use of the implant,

INDICATION OF USE - SILIMED Silicone Gel Mammary Implants are indicated for women for the following procedures:

- Breast enlargement. Breast augmentation includes preliminary breast augmentation to increase breast size, as well as revision surgery to correct or improve the outcome of preliminary breast augmentation surgery.

- Breast reconstruction. Breast reconstruction includes primary reconstruction to replace breast tissue that was removed due to cancer or trauma or that failed to develop properly due to a serious breast abnormality. Breast reconstruction also includes revision surgery to correct or improve the results of a primary breast reconstruction surgery.

CONTRAINDICATIONS - Groups of patients in which the implantation procedure is contraindicated:

- Patients with infection anywhere in the body;
- Patients with cancer or pre-malignant condition not treated properly;
- Pregnant or breastfeeding patients;
- Patients with atopy;
- Patients with suppressed immune system;

¹ High Strength Silicone

- Patients with local inflammation.

PRECAUTIONS - The physician must, in the preoperative evaluation of his patients, take into account: recent history of tumors in the region of the implant, in particular recurrent cancer or metastasis; autoimmune diseases; pathologies and treatments that affect the immune system or blood clotting and; history of severe allergy; diabetes and cardiovascular problems.

INSTRUCTIONS FOR HANDLING AND USE -

- Package Opening:

1 - Make sure that the outer plastic cover has not been opened;

2 - Remove the double packaging from inside the sealed box. It must be carefully examined before use in the operating room. The product that has its packaging violated in any way must not be used;

3 - Separate the documents that come with the product;

4 - Stick the adhesive labels with the product data to the hospital, doctor's files and also to the dossier that will be delivered to the patient, as indicated by the symbols;

5 - Open the outer blister to have access to the sterile inner blister that contains the product. Be careful not to contaminate it outside the first one;

6 - Open, then, the sterile internal blister, in the surgical field.

The high dielectric value of silicone can generate electrostatic charges, which are responsible for attracting particles in the environment, such as fluff in general and talc, to mention a few examples. External contaminants that adhere to the shell surface can generate reactions such as those caused by foreign bodies to the body, through the increase of fibrosis and the generation of fluids. For this reason, it is very important to be careful to open the packaging.

SURGICAL PROCEDURE - Adequate surgical techniques and procedures are necessarily the responsibility of the physician, who must evaluate them according to his own training and experience, observing updated and accepted techniques to minimize the occurrence of adverse reactions. The most common techniques and complications are described in this document. In addition, SILIMED has a scientific bibliography that can be requested through its representatives or directly to the headquarters.

It is also the physician's responsibility to make sure that his team is properly trained in order to perform the required medical procedures.

The physician should discuss the potential risks and benefits of a preoperative mammography and ultrasound for future reference. In case of suspicious changes in these exams, we suggest the investigation before or during the insertion of the implant.

Before the surgical procedure, the patient must be informed of all the procedures to which she will be submitted and all the risks involved. The application of preoperative informed consent is recommended.

The physician must carefully evaluate the type of implant, access route, dissection of the pocket and location of the implant, taking into account the patient's anatomy and the clinical and aesthetic requirements of each case.

Because they are filled with high-cohesive silicone gel, SILIMED implants can be introduced through 3 access incisions: inframammary, periareolar and transaxillary. The choice of incision should be in accordance with the patient's anatomy and the surgeon's experience.

The inframammary route is technically the simplest, due to the easy access to the subglandular, subfascial and submuscular planes. It is also easy to make the pocket, perform hemostasis, insert the implant and confirm its positioning.

The half-moon,shaped periareolar access demands a dissection bordering the limit of the mammary gland or sectioning it until reaching the subglandular, subfascial and submuscular planes. Depending on the thickness of the gland and the surgeon's experience, it may be simpler or more laborious to make the pocket in the chosen plane, perform hemostasis and introduce and confirm the implant positioning. The creation of an asymmetric pocket between the breasts can lead to an asymmetry of the breast groove. Periareolar introduction may also be associated with a higher incidence of infection and biofilm formation due to the contact of the implant with the ruptured breast ducts during the introduction of the implant, in addition to increasing the risk of interference with breastfeeding.

Transaxillary access requires the dissection of a tunnel between the axillary incision and the mammary space. Therefore, the distance between the axillary access and the limits of the pocket regardless of the plane - subglandular, subfascial and submuscular - is greater. Some surgeons may find technical difficulty in performing hemostasis and precise dissection of the limits of the mammary sulcus by this incision. The introduction of the implant requires greater skill and care in its introduction, in order to avoid trauma to the axillary skin and fracture of the implant gel due to excessive force. The incorrect positioning of the implant can lead to groove asymmetry, rotation and displacement of the implant.

Regardless of the incision chosen, the skin incision must be appropriate for the correct visualization and dissection of the pocket. We suggest that the incision be between 5 and 7 cm, depending on the type of access chosen and the size of the implant. The disproportion between these aspects can lead to excessive force causing burns at the edges of the incision due to friction or fracture of the implant gel. (Excessive pressure) should be avoided during implantation.

With regard to positioning, three allocations for the implant can be used: retroglandular, retrofascial and retromuscular. All have advantages and disadvantages. These must be exposed and explained to the patient before implantation. As a general rule, the pre-pectoral muscle pathways (retroglandular and retrofascial) preserve the pectoral musculature and have a lower incidence of postoperative pain. However, in patients submitted to the inclusion of large implants or small breasts, the implants, especially their edges, can be palpable and visible. The retromuscular plane, on the other hand, has a thicker coverage, making palpation and visualization of the implant difficult, but it is more painful, changes muscle function and can lead to poor positioning of the implant with lateral and superior displacement.

Once the implants are introduced, one must confirm their

position. In round and tapered implants, advance, maximum models, the positioning symmetry between the breasts must be determined, avoiding poor positioning in the postoperative period, with lateral, cranial or caudal displacement of the implant, as well as groove asymmetries. In the natural, nuance and enhance model implants, in addition to confirming the symmetry of position between the breasts, as previously described, make use of the colored implant marking (visual indication) and high relief markings on the implant base (tactile indication), to determine its correct orientation, avoiding positions with rotated or asymmetrical axes.

The plastic glove accessory aims to assist the inclusion of the mammary implant - textured surface in the implantation pocket. The implant must be placed inside the device through one of the openings and, in the doctor's manipulation, the implant is introduced into the pocket. After ensuring the placement of the implant in the desired position, this accessory must be removed.

It is known that the greatest immediate perioperative risks are: hematoma, seroma and infection. Therefore, correct asepsis, hemostasis and implant handling are recommended in order to reduce these complications. The use of drains is possible according to the surgeon's preference and / or the patient's need.

Basic postoperative recommendations include wearing the appropriate bra for the procedure performed, avoiding activities that may favor pain, bleeding or displacement of the implants, and avoiding smoking.

For postoperative follow-up, the evaluation of the patient and the implant must be done through anamnesis and physical examination, in order to detect possible changes in the implant. Imaging tests such as ultrasound and nuclear magnetic resonance are the best to date to assess the integrity of the implants. We suggest performing an annual ultrasound to evaluate the implants.

The physician should always have a backup product available during surgical procedures.

We suggest an annual clinical follow-up.

INSTRUCTIONS FOR REMOVAL OF THE IMPLANT - The removal of mammary implants, in general, represents a simple procedure. However, possible previous complications can add some complexity to the surgery.

The access incision for the removal of the implant must be the most appropriate at the physician's discretion, based on previous exams. After incision in the skin and dissection of the tissue planes, the plan used to place the implant must be accessed. Proceed with careful removal of the implant. Then close the layers and finish the surgical procedure.

It is up to the doctor to choose whether or not to remove the fibrotic capsule, assessing the impact of this removal on the patient's morbidity in the postoperative period and also on possible interference in the interpretation of mammographic exams. In case of permanent removal, without subsequent replacement of the implant, the doctor must also evaluate the aesthetic result that the extraction of the capsule will cause.

In case of removal of the implant whose membrane has

suffered loss of integrity due to contact of the gel with the patient's tissues, the doctor must pay attention to the removal of the entire gel content, and must also proceed with the removal of the fibrotic capsule.

Meticulous hemostasis during surgery is recommended in order to decrease the chances of postoperative hematoma and seroma. Excessive and persistent bleeding must be controlled throughout the procedure.

The damaged product must be returned to the manufacturer for disposal.

IMPORTANT RECOMMENDATIONS / WARNINGS -

Mammary Implants can only be purchased by doctors or under prescription.

It is supplied sterile and, to maintain its aseptic condition during surgery, asepsis and cleaning in the conditions of its use is necessary.

The implant and its packaging must be intact, otherwise the implant must not be used. It can only be used if found intact in its original form of manufacture, that is, without any change in its original characteristics.

Eventually, the implant may show bubbles or a certain whiteness, resulting from sterilization, which does not prevent its use.

The implant surface must not be contaminated with talc, dust or oils. Before contacting the implant, the medical team must replace the sterile surgical gloves in use with new pairs, and wash them immediately with sterile and apyrogenic saline.

Solutions containing iodine cannot come into contact with the implant.

Marker solutions, such as methylene blue, used in the medical-surgical procedure must be properly sterilized to avoid biological contamination.

All perforating surgical instruments to be used in the region close to the implant must be handled with care, so that it is not affected. Any puncture, cut or accidental scratch on the implant shell will invalidate it for surgical use. Damaged product repair is not possible. It must be replaced with another one.

The physician should instruct his patient not to massage and / or vigorously manipulate the implant or to engage in intense physical activity for a specified period.

The patient should be instructed to follow the most up-to-date medical recommendations regarding breast exams and to distinguish the implant from the breast tissue, in order to obtain the best result with self-examination.

Although there is no definitive data yet, the scientific literature suggests that there is no direct correlation of cause and effect between the placement of breast implants for aesthetic increase and suicide frequency, and the increase in this frequency is possibly correlated with intrinsic psychiatric morbidity of these patients. Therefore, it is highly recommended that the psychological stability of the patient is carefully evaluated by the doctor before and after the placement of Mammary Implants.

1- Breast implants are not considered lifelong devices. The longer people have them, the more likely they are to develop complications, some of which will require more surgery.

2 - Breast implants have been associated with the development of a cancer of the immune system called anaplastic large cell lymphoma (BIA-ALCL). This cancer occurs more commonly in patients with textured breast implants than smooth implants, although the rates are not well defined. Some patients died of BIA-ALCL.

3 - Patients who received breast implants reported a variety of systemic symptoms such as muscle and joint pain, confusion, chronic fatigue, autoimmune diseases, among others. The patient's individual risk of developing these symptoms has not been well established. Some patients report complete resolution of symptoms when these implants are removed without replacement.

The Mammary Implant is intended for single use. Prohibited reprocessing.

By law, it cannot be reused or resterilized, as it damages the product and compromises its performance and safety.

The damaged product must be returned to the manufacturer for disposal.

POTENTIAL COMPLICATIONS AND POSSIBLE SOLUTIONS

- Any patient who has undergone surgery to introduce a foreign element to the body is subject to possible complications.

SILIMED delegates to physicians the responsibility to inform patients about the possibility of new surgery for removal or replacement of the implant, as well as the possible occurrence of adverse effects after its implantation and the appropriate treatment. At the same time, it invites the medical profession to communicate other findings in this regard.

SILIMED also delegates to the doctor to adopt the best solutions available for the complications of surgery. Among these, it is worth mentioning the meticulous hemostasis during the procedure, antibiotic therapy and removal or replacement of the implant.

COMPLICATIONS PREDICTED - For the Mammary Implant, the potential complications of the surgery, reported in the current medical-scientific literature, or directly to SILIMED, followed through the product risk management process - post-production phase, are as follows:

CHEST-WALL FLATTENING- Few cases have been reported related to tissue expanders and, occasionally, to mammary implants placed in a retromuscular position. It is associated with the surgical technique used. It can cause pain and injury to the organs near the ribs and muscles of the thoracic region.

ALLERGENICITY (SKIN ERUPTION - RASH) - It can occur in the recent postoperative period and is mainly characterized by tissue irritation, redness and rash. It is easily treatable and of temporary duration.

CHANGE IN BREAST AND / OR NIPPLE SENSITIVITY - Uncommon. Corresponds to an increase (hyperesthesia) or decrease (hypoesthesia) in the sensitivity of the implanted

region due to the surgical incision. This is usually a process that is reversible over time, but it can take two years or more.

TISSUE ATROPHY - Associated with large volume mammary implants, when inserted in pockets disproportionate to the volume of these implants. Mechanical forces press on the parenchyma and subcutaneous tissue, which can cause atrophy of these tissues over time.

CALCIFICATION - It is characterized by the presence of deposits of insoluble calcium (and / or magnesium) salts, especially calcium carbonate and phosphate in the fibrotic capsule. It can cause pain and hardening of the scar tissue. Of unknown origin, it was observed in the tissues surrounding the implantation, sometimes requiring removal of the implant.

CYSTS / NODULES - Breast nodules and cysts are changes that happen in breast tissue. Breast lumps are breast protuberances. They can usually be caused by fibrocystic changes, fibroadenomas, cysts and cancer. They can also occur due to rupture of the implant, if leakage occurs and the silicone gel migrates to lymph nodes or other parts of the body. The cyst is characterized as a closed sac / cavity filled with liquid covered by an epithelium.

IATROGENIC COMPLICATIONS - Certain situations such as folds and perforations in the mammary implant, among others, can be caused by the use of inappropriate surgical technique.

CAPSULAR CONTRACTURE / FIBROSIS - Occurs when the scar tissue or fibrotic capsule becomes thicker and retractable and compresses the implant. The fibrotic capsule forms around the implant due to the body's reaction to the presence of a foreign body. Factors such as infection, hematoma and seroma can favor the appearance of capsular contracture. The probability of occurrence decreases with meticulous hemostasis of the pocket and care in handling the implant. SILIMED advises against the external (closed) capsulotomy procedure, which can cause the implant to rupture.

SUTURE DEISCENCE - Suture dehiscence is the separation of the layers of the surgical wound. It can be partial and superficial or complete, when there is complete rupture of all layers of the tissue. This effect usually starts with a small seroma. In most cases, this is a problem related to the surgical technique and / or the use of an implant that is larger than adequate.

DISPLACEMENT OF THE IMPLANT / BAD POSITIONING - It can happen with patient discomfort and distortion in the desirable way, this being a problem related to the operative technique, ineffective formation of the fibrotic capsule or the occurrence of trauma. The physician should assess the need for removal of the implant and / or its likely replacement

PAIN - Some patients may experience pain in the postoperative period, which should be investigated immediately. It can occur from mild to severe and from short to long lasting. It can result from infection, inflammation, seroma, calcification and / or the use of inappropriate surgical technique.

WRINKLING - There may be a ripple on the surface of the mammary implant that can be seen and / or felt.

EROSION - After implantation, tissue erosion near the implant can occur, caused by folds in the device and other deformities that lead to exaggerated compression of the tissue.

EXTRAVASATION OF THE SILICONE GEL - This is the process of diffusion and migration of the gel by rupture of the implant to the adjacent tissues. It can cause the appearance

of complications, such as: granuloma, allergy, cysts and / or nodules, immune response, pain, edema, breast hardening, flaccidity, tissue infiltration, irritation, capsular contracture, swelling of the lymph nodes. The SILIMED high performance silicone gel minimizes the occurrence of this adverse effect.

EXTRUSION - Unstable or compromised tissue covering, excessive stress or trauma occurring at the implantation site and / or interrupting the healing process may result in the extrusion of the implant. It has been reported as an uncommon complication in both patients who have undergone augmentation and reconstruction surgery.

GALACTORREA - May be the result of hormonal disorders or the action of medications. It causes discomfort and embarrassment, affecting the patient's emotional state.

GRANULOMA - It is associated with the formation of granular tissue due to the presence of several types of cells such as macrophages, lymphocytes and fibroblasts. It can lead to a chronic inflammatory lesion, tumor and consequent surgical intervention to remove the implant.

HEMATOMA - The extravasation of blood and its confinement in an organ or tissue (hematoma) can occur in any type of surgery. In its occurrence, it must be carefully evaluated to define the conduct to be adopted. With the risk of being a precursor to infection and increased fibrosis, if it is not reabsorbed, the hematoma must be removed. The main preventive measure is perfect hemostasis during surgery and postoperative rest.

INFECTION (BIOLOGICAL CONTAMINATION) - Uncommon and difficult due to the implant, if the recommendations for use and asepsis associated with implantation surgery are observed. It includes the possibility of secondary necrosis and extrusion of the implant. Local culture is recommended and, in sequence, the use of antibiotics, drainage and, eventually, removal of the implant.

INTERFERENCE WITH MAMMOGRAPHY, ULTRASOUND AND MAGNETIC RESONANCE - Interference occurs when the presence of implants (or artifacts) causes an erroneous reading of the result, for example: false positive diagnosis of rupture, resulting in unnecessary removals. Mammary implants can interfere with the diagnosis of breast cancer made by mammography, covering up any suspicious spots on the breast tissue, masking the results and delaying the start of anti-cancer treatment.

INTERFERENCE IN BREASTFEEDING - The presence of silicone in breast milk has been researched in several studies, with no significant levels of the substance found in mothers with silicone mammary implants, when compared to women without implants. However, there is a possible interference in breastfeeding (reduced milk production) in women who have undergone mastopasty with mammary implants. The frequency of interference with breastfeeding is higher in patients in whom the breast implant incision was periareolar.

LARGE CELL ANAPLASIC LYMPHOMA ASSOCIATED WITH BREAST IMPLANT (BIA-ALCL) - BIA-ALCL is a type of non-Hodgkin's lymphoma, which is a cancer of the immune system (not breast cancer), most often found in tissue scar that surrounds the implant, also called a capsule, or in the fluid that forms between the implant and the capsule (seroma). The occurrence of late seromas is widely considered by the plastic

surgery community to be associated with the pathogenesis of ALCL. The pathogenesis of BIA-ALCL is still unknown and its development is slow. The mechanism believed to be related to the development of BIA-ALCL is the body's inflammatory response (foreign body reaction) in an exacerbated way. Another pathway that can trigger BIA-ALCL is bacterial infection by biofilm around mammary implants, leading to an inflammatory response. An important factor that may be associated with the development of the disease is the genetic issue. Although an increase in the prevalence of BIA-ALCL has been observed, its incidence is currently low. However, although any type of mammary implant, regardless of surface, filling or manufacturer, can be associated with the development of ALCL, the prevalence is higher in implants with a textured surface.

Different treatments have been reported, such as implant removal and capsulotomy, lymph node dissection, chemotherapy, radiotherapy and even autologous stem cell transplantation. However, there is no universal treatment. If diagnosed, individual treatment should be performed. Women should monitor their mammary implants and contact their doctors if they notice any changes, and those considering mammary implant surgery should discuss the risks and benefits with their doctors.

NECROSIS - Resulting from exaggerated skin tension caused by the presence of the implant or by surgical trauma. The physician should assess the need for removal of the implant and / or its likely replacement.

LOSS OF MECHANICAL INTEGRITY / PRODUCT FAILURE - In case of loss of mechanical integrity of the mammary implant, the patient may be submitted to surgery for replacement, depending on the degree of danger. Among these situations, we can mention the detachment of the implant pellet, the spontaneous rupture of its shell, among others.

PNEUMOTHORAX - This is the collapse of the lung due to a sudden change in pressure within the chest cavity. In the case of mammary implants, this can occur when the implantation pocket is formed due to surgical trauma, especially in the retromuscular space. The consequences are: tachycardia, chest pain, cough, shortness of breath, cyanosis, and agitation, among others.

INFLAMMATORY RESPONSE - The presence of an element foreign to the body triggers an inflammatory process, which can be more or less intense. Intense inflammation can also be due to trauma or inadequate procedures. It causes pain, redness, local heat and swelling. This effect can be treated pharmacological

INSATISFACTORY AESTHETIC RESULT - It is associated with the patient and / or doctor's dissatisfaction with the result of the mammary implant placement surgery. Some possible unsatisfactory aesthetic results are: ptosis, flaccidity, asymmetry, folds, folds, implant wrinkling, stretch marks, adhesion and hypertrophic scars. The doctor must evaluate the best measure to be taken, without necessarily having to undergo a surgical intervention to change the implant.

RUPTURE - It can be generated by the natural wear of the implant due to the implantation time, by damage by a sharp instrument, by excessive effort during surgery, by trauma or extreme physical pressure, by compression during mammography or by severe capsular contracture. Most of them occur during the surgical procedure. In case of suspected rupture of the implant in an asymptomatic patient ("silent

rupture”), a careful assessment should be made to confirm the situation and decide on the procedures to be taken, including the possibility of replacing the implant.

SEROMA - Accompanied by edema and pain, which may subside with rest, immobilization and application of compresses, eventually drainage and even removal of the implant. The seroma is not always acute. There may be a minimal, continuous, undetected seroma, which appears until months later as “late seroma”. There are still “late seromas” of unknown origin.

SYMPTOMS ASSOCIATED WITH MACROMASTY - They are associated with large volume breasts. Some of the symptoms are: back pain, breasts, hands (including numbness), neck, head; skin irritation and shoulder marks. The volume of the mammary implant must be considered respecting the amount of breast tissue present, the shape of the breasts and the patient’s desire. The methodologies used to estimate whether the volume of mammary implants is adequate to the needs and expectations of patients are of great importance when choosing a mammary implant.

GEL TRANSUDATION (GEL BLEEDING) - It is associated with mammary implants filled with silicone gel. It is the process of diffusion and migration of the gel through the intact membrane (transudation) or by the rupture of the implant to the adjacent tissues. It can cause the appearance of local complications, such as: reduction in size, asymmetry, change in the shape of the sinus, pain, capsular contracture or formation of nodules and granulomas. The SILIMED high performance silicone gel minimizes the occurrence of this adverse effect.

THROMBOSIS - Uncommon, usually associated with surgical procedure or with patients with obesity, neoplasms, varicose veins, traumatic or iatrogenic venous lesions, patients bedridden or immobilized by a systemic or local disease in the limb, patients with heart disease, extrinsic vein compressions, pregnancy and puerperium, use of hormonal contraceptives and hematological diseases. It can cause the appearance of local complications, requiring the removal of the implant.

Other complications: Emotional change, Discomfort, Numbness, Edema, Erythema, Breast hardening, Hemorrhage, Swollen lymph nodes (axillary lymph nodes) also called Lymphadenopathy, Tissue infiltration, Irritation, Tissue injury, Palpability, Itching, Flushing, Skin tension, Implant visibility, Asymmetry, Ptosis, Late healing, Upper pole filling.

Some complications have been reported in the medical-scientific literature, but so far there is no scientific evidence of a possible cause-and-effect relationship between these events and silicone implants. Some examples are: **BREAST CANCER** - According to current medical-scientific literature, women undergoing silicone implant surgery are not at a higher risk of developing breast cancer.

IMMUNOLOGICAL RESPONSE - Studies carried out to date have found no evidence that the placement of silicone implants causes connective or autoimmune diseases. The current literature suggests a possible link between implant rupture and fibromyalgia and connective tissue diseases. Further studies are needed.

BREAST IMPLANT ILLNESS - Recent studies have focused attention on a set of systemic symptoms, which are known as “Breast Implant Illness”. This term appeared mainly on social media, and represents an extensive group of poorly defined

symptoms, which would be related to breast implants. Over the years, some syndromes, such as human adjuvant disease (1964), silicone-induced human adjuvant disease, adjuvant-induced autoimmune / inflammatory syndrome (2011) and silicone implant incompatibility syndrome (2013) have been used to link diseases systemic to silicone and other adjuvants. In general, the term BII was applied loosely to refer to these entities and, sometimes, more comprehensively, to include all complications related to silicone breast implants. Due to the lack of scientific evidence and the lack of evidence of the correlation of these symptoms with the presence of the implant, breast implant disease is not yet recognized as a disease itself.

TERATOGENICITY OR OTHER REPRODUCTIVE EFFECTS - It is the ability to cause congenital malformations in the fetus. The review of the medical-scientific literature indicates that studies have shown no evidence of teratogenicity or other reproductive effects associated with mammary implants.

EFFECTS ON CHILDREN BORN TO MOTHERS WITH MAMMARY IMPLANTS - Recent studies have shown that the amount of silicone in milk found in women who have silicone implants is not significantly different from women who do not have silicone implants. However, the review of the medical-scientific literature does not demonstrate evidence that silicone gel breast implants can cause any harmful effects to the children of women with implants.

SUICIDE - Epidemiological studies indicate that women who have undergone cosmetic augmentation mammoplasty have a higher suicide rate than the general population. Although there is still no definitive data, the scientific literature suggests that there is no direct correlation of cause and effect between the placement of mammary implants for aesthetic increase and suicide frequency, the increase in this frequency being possibly correlated with the intrinsic psychiatric morbidity of the patients. Therefore, it is highly recommended that the psychological stability of the patient is carefully evaluated by the doctor before and after the placement of mammary implants.

DURABILITY - Like all breast implants, the SILIMED Mammary Implants have a limited useful life. This lifetime has not been determined by the scientific community yet. SILIMED establishes a mean period of 10 years as the expected useful life, but such parameter may be altered whenever some reason to justify it is assessed.

The implant may need to be removed or replaced, which may imply in a revision surgery.

An annual clinical follow-up is recommended for these patients.

INFORMATION FOR SILIMED PATIENT - The following information is intended for all patients who intend to undergo augmentation or reconstructive surgery using silicone mammary implants. They are information about the main doubts about the subject, including potential complications, benefits and risks, activities that can damage the implant, the possible need for removal / replacement of devices and the expected durability for the mammary implant.

The physician responsible for the procedure should carefully read this section, being responsible for resolving all the patient’s doubts before the procedure with SILIMED mammary implants.

Annual clinical follow-up of these patients is recommended.

If you suspect any complications listed in this newsletter or any abnormal symptoms, contact your doctor.

EXPECTED BENEFITS - Mammary implant surgery can be of great benefit in providing successful breast reconstruction, augmentation and surgical revision. Reconstruction / augmentation with a mammary implant also results in improvements in quality of life. In addition, studies show that breast reconstruction with a mammary implant has been an aid in the recovery of breast cancer, in addition to reducing emotional stress by helping the body to return to a more natural appearance, as opposed to not having reconstructive surgery, or use external prosthesis.

CLARIFICATION AND CONSENT OF THE PATIENT -

Considering the risks inherent to a surgery, with or without the use of implants, and the possible related consequences, SILIMED relies on the physicians to clarify their patients as to the existing risk-benefit balance.

Every patient shall be given the SILIMED Information to the Patient during the initial consultation so that the patient is allowed sufficient time prior to the surgery to suitably read and understand the important risk information, follow-up recommendations and benefits associated to the surgeries for placement of mammary implants. For a successful documentation of the patient's clarification, the Clarification Term for use of the SILIMED Mammary Implant shall be signed by the patient and by the physician and then attached to the patient's medical register.

INTERFERENCE WITH MAMMOGRAPHY; UNTRASOUND AND MAGNETIC RESSONANCE - Interference occurs when the presence of implants (or artifacts) causes an erroneous reading of the result, for example: false positive diagnosis of rupture, resulting in unnecessary withdrawals.

SILIMED points out that the implant can interfere with the quality of the mammogram. For this reason, the patient must be guided and instructed to request the services of professionals who have experience with mammography techniques with implant, informing the radiologist of the need to adapt the mammographic compression and not forgetting to show the "PRODUCT TRACEABILITY CARD SILIMED - FOR PATIENT".

Other methods, such as ultrasound and magnetic resonance imaging, may be useful in conjunction with mammography because they do not require compression and allow examination from any angle.

According to recent studies, silicone breast implants reduce the image quality of echocardiographic visualizations. More data is needed, however SILIMED advises its patients to inform the technician responsible for the exam about the existence of silicone implants.

CAUTION - In the first month after surgery, some activities can damage the implant and should be avoided, such as:

- Exposure to the sun;
- Sudden movements;
- Sports in general.

You should be aware that normal stress or trauma applied to the surgery site can cause the extrusion of the implant.

Wearing an appropriate bra and not doing violent exercises are minimum recommendations to be followed in the postoperative period.

The patient must inform the medical or pharmacist specialist that she has silicone implants before using topical medication in the breast region.

The patient needs to continue to consult with the medical specialist to proceed with routine monitoring to detect breast cancer.

The patient should inform the doctor about the presence of the implant if any breast surgery is scheduled.

If you suspect any complications listed in this newsletter or any abnormal symptoms, especially if trauma or compression occurs (e.g. excessive breast massage through some sports or wearing a seat belt), contact your doctor.

IMPORTANT FACTORS TO CONSIDER -

1. If you will undergo breast augmentation or reconstruction, be aware that breast implants are not considered devices for life and that breast implantation may not be a one-time surgery. It is possible that you may need one or more additional surgeries and have to schedule visits to the doctor throughout your life. It is also possible that you may need surgery to remove the implant, with or without replacement.

2. Many of the changes in your breast after implantation cannot be undone. If you choose to have your implant(s) removed in the future, you may experience unacceptable folds, drooping or wrinkling, loss of breast tissue or other alterations to your breast contour.

3. With breast implants, routine screening mammography will be more difficult and you will need additional images, which means more radiation and more time.

4. Every SILIMED Breast Implant comes with a "Silimed Product Traceability Card - For Patient", which contains all information regarding the surgery. This card is for your safety, you must carry this card to facilitate medical care in an emergency.

STORAGE / CONSERVATION / TRANSPORTATION - SILIMED products and packages are resistant and if stored at environment temperature and according to the guidance printed in the package, they will not be damaged. Since they are sterile medical products, it is recommended that they are not transported or stored along with other types of materials which can cause physical or chemical damages to the product, which would prevent its use.

WARRANTY - SILIMED will replace any of its line products that presents a provable manufacturing defect, provided it is returned duly identified by the person who purchased it and after analysis at the SILIMED laboratory.

The adverse effects presented here, the inappropriate surgical use and the color variation mentioned in the item "Constituent Materials", are not considered product defects.

The guarantee for SILIMED products does not cover the simple decision of the patient or the doctor to change the product.

TRACEABILITY - SILIMED maintains an accurate register, unit

by unit, of raw materials, manufacturing stages, environmental and operational conditions and quality tests, assigning a single number to each product that completely identifies it at any time. In this sense, it is absolutely necessary that any complaint shall be sent together with the product's SERIAL NUMBER (SN) or, when applicable, LOT NUMBER, of the unity complained, contained on the packaging.

SILIMED provides adhesive labels with the data on the product, which shall be attached to the hospital's card and to the physician's card, as indicated by the symbols.

Also provided are the Patient's Registration Card, to be signed by the patient and returned to SILIMED by the physician; and the Card for Traceability of SILIMED Product, to be handed to the patient by the physician. These documents contain the data on the products, which allow their traceability.

RETURN OF PRODUCTS - Products contaminated with biological material shall be properly decontaminated by the physician or responsible health unit before being sent back to SILIMED, in accordance with the applicable legislation.

No product shall be returned without the prior authorization by SILIMED.

Dear Doctor: In order to facilitate understanding when handling our products, in this printed notice you will find international symbols, together with those developed at SILIMED, that appear on our labels, packaging, instructions for use and/or catalogues, as well as their meanings.

	Single use		Do not use if package is damaged
	Lot number		Volume of the product
	Date of manufacture		Consult instructions for use
	Ethyleneoxide sterilization		Keep away from sunlight
	Steam or dry-heat sterilization		Dimensions of product
	Manufacturer		CE mark
	Use by		Quantity of product
	Serial number		ANVISA registration number
	Sterile		Label for tracing to the doctor's dossier
	Reference number		Verified by the QC

	Caution, consult accompanying documents		INMETRO seal for identification of conformity
	Authorised representative in the European community		Label for tracing to the hospital dossier
	Label for tracing the patient's dossier		Medical device
	Keep dry		Temperature limit
	This way up		Humidity limit
	Right Side		Unique Device Identifier
	Left Side		Portuguese
	Spanish		Italian
	Russian		English

EN

INFORMATION FOR PATIENT - SILIMED MAMMARY IMPLANTS

Silimed has developed this leaflet to help you understand the benefits and risks related to augmentation or reconstructive surgery using silicone breast implants. After you have read this leaflet carefully, talk to your doctor before making any decisions. If you decide on the procedure, you and your doctor must sign the Informed Consent Form prior to surgery, confirming that you have read and understood the information provided about the benefits and risks of surgery with silicone breast implants. The Clarification Term for Use of the SILIMED Mammary Implant is at the end of this leaflet.

DESCRIPTION - SILIMED Mammary Implants consist of a single outer elastomeric shell, filled with medical grade silicone gel; they can have a smooth, textured surface or even coated with polyurethane foam. They are all supplied sterile and intended for single use.

These implants are available in different sizes and shapes to help each woman achieve the best result for her body. The doctor should guide you in choosing the correct implant to obtain the desired aesthetic effect.

EXPECTED BENEFITS – Breast augmentation / reconstruction surgery can bring great benefit by providing reconstruction, enlargement and surgical revision of the breast with a high level of satisfaction. Patients choose primary breast augmentation surgery to increase the size and proportion of their breast (s). In addition, patients choose revision surgery (replacement of an existing breast implant) to correct or improve the results of primary augmentation surgery. According to the medical-scientific literature, most women who have undergone surgery with breast implants have reported high levels of satisfaction with body image, shape and breast size. Additionally, after surgery, patients reported greater self-esteem, and felt more feminine and attractive. Reconstructive surgery with a breast implant also results in improved quality of life. In addition, studies demonstrate that breast reconstruction with a breast implant has been an aid in the recovery of breast cancer, in addition to reducing emotional stress by helping the body to return to a more natural appearance, as opposed to not having surgery reconstructive or use external prosthesis.

CLARIFICATION AND CONSENT OF THE PATIENT - Considering the risks inherent in surgery, with or without the use of implants, and the possible complications related to them, SILIMED counts on doctors to clarify their patients regarding the existing risk-benefit balance.

Each patient should receive the Patient Information - SILIMED Mammary Implants during their initial consultation to allow sufficient time before surgery to properly read and understand important risk information, follow-up recommendations and benefits associated with breast implant placement surgeries. To successfully document the patient's clarification, the SILIMED Mammary Implant Clarification Term must be signed by the patient and the doctor and attached to the patient's medical record.

EXPECTED COMPLICATIONS- Breast implant surgery offers great benefits. However, any type of surgery involves risks. There are some local complications (problems near the surgery site) that can occur after surgery. SILIMED delegates to physicians the responsibility to inform patients about the possible occurrence of adverse effects, consequences and appropriate treatment. Below, we present the potential risks and complications, reported in the literature, associated with surgery with breast implants of different models and manufacturers, as well as the possible effects related to these risks.

Alteration of Breast and / or Nipple Sensitivity, Chest-wall flattening, Allergenicity / Rash - Rash, Emotional alteration / Depression, Tissue Atrophy, Calcification, Cysts / nodules / tumor, Iatrogenic complications, Capsular contracture / fibrosis, Suture dehiscence , Discomfort, Implant Dislocation, Pain, Numbness, Edema (Swelling), Erythema (Hyperemia), Breast Hardening, Wrinkling, Erosion, Extrusion, Galactorrhea, Granuloma, Hematoma, Hemorrhage, Implant Removal, Infection (Biological Contamination) and its symptoms, swollen lymph nodes (axillary lymph nodes), tissue infiltration, irritation, tissue injury, interference with breastfeeding, poor positioning, silicone migration (lymphadenopathy and silicoma), necrosis, loss of mechanical integrity / product failure, palpability, pneumothorax, Pruritus, Local Inflammatory Response, Revision Surgery, Unsatisfactory Aesthetic Result (such as ptosis, flaccidity, asymmetry, folds, adhesion , stretch marks and hypertrophic scars), Flushing, Rupture (including cases of "silent" rupture), Secretion, Seroma / fluid, Symptom, Symptoms Associated with Macromastia, Skin tension, Gel transudation (gel bleeding), extravasation of the silicone gel, thrombosis, vesicles, implant visibility, and filling of the upper pole.

LARGE CELL ANAPLASIC LYMPHOMA ASSOCIATED WITH BREAST IMPLANT (BIA-ALCL) - BIA-ALCL is a type of non-Hodgkin's lymphoma, which is a cancer of the immune system (not breast cancer), most often found in tissue scar that surrounds the implant, also called a capsule, or in the fluid that forms between the implant and the capsule (seroma). The occurrence of late seromas is widely considered by the plastic surgery community to be associated with the pathogenesis of ALCL. The pathogenesis of BIA-ALCL is still unknown and its development is slow. The mechanism believed to be related to the development of BIA-ALCL is the body's inflammatory response (foreign body reaction) in an exacerbated way. Another pathway that can trigger BIA-ALCL is bacterial infection by biofilm around mammary implants, leading to an inflammatory response. An important factor that may be associated with the development of the disease is the genetic issue. Although an increase in the prevalence of BIA-ALCL has been observed, its incidence is currently low. However, although any type of breast implant, regardless of surface, filling or manufacturer, can be associated with the development of ALCL, the prevalence is higher in implants with a textured surface.

Different treatments have been reported, such as implant removal and capsulotomy, lymph node dissection, chemotherapy, radiotherapy and even autologous stem cell transplantation. However, there is no universal treatment. If diagnosed, individual treatment should be performed. Women should monitor their breast implants and contact their doctors if they notice any changes,

and those considering mammary implant surgery should discuss the risks and benefits with their doctors.

Other complications have been reported in the medical-scientific literature, however, to date, there is no scientific proof of a possible cause and effect relationship between these events and silicone implants. Some examples are:

BREAST CANCER - According to current medical-scientific literature, women undergoing silicone implant surgery are not at a higher risk of developing breast cancer.

IMMUNOLOGICAL RESPONSE - Study To date, no evidence has been found that the placement of silicone implants causes connective or autoimmune diseases. The current literature suggests a possible link between implant rupture and fibromyalgia and connective tissue diseases. Further studies are needed.

BREAST IMPLANT ILLNESS - Recent studies have focused attention on a set of systemic symptoms, which are known as “Breast Implant Disease”. This term appeared mainly on social media, and represents an extensive group of poorly defined symptoms, which would be related to breast implants. Over the years, some syndromes, such as human adjuvant disease (1964), silicone-induced human adjuvant disease, adjuvant-induced autoimmune / inflammatory syndrome (2011) and silicone implant incompatibility syndrome (2013) have been used to link diseases systemic to silicone and other adjuvants. In general, the term BI was applied loosely to refer to these entities and, sometimes, more comprehensively, to include all complications related to silicone breast implants. Due to the lack of scientific evidence and the lack of evidence of the correlation of these symptoms with the presence of the implant, breast implant disease is not yet recognized as a disease itself.

TERATOGENICITY OR OTHER REPRODUCTIVE EFFECTS - It is the ability to cause congenital malformations in the fetus. The review of the medical-scientific literature indicates that studies have shown no evidence of teratogenicity or other reproductive effects associated with breast implants.

EFFECTS ON CHILDREN BORN TO MOTHERS WITH BREAST IMPLANTS - Recent studies have shown that the amount of silicone in milk found in women who have silicone implants is not significantly different from women who do not have silicone implants. However, the review of the medical-scientific literature does not show evidence that silicone gel breast implants can cause any harmful effects to the children of women with implants.

SUICIDE - Epidemiological studies indicate that women who have undergone cosmetic augmentation mammoplasty have a higher suicide rate than the general population. Although there is still no definitive data, the scientific literature route suggests that there is no direct correlation of cause and effect between the placement of breast implants for aesthetic increase and suicide frequency, the increase in this frequency being possibly correlated with their intrinsic psychiatric morbidity. patients. Therefore, it is highly recommended that the psychological stability of the patient is carefully evaluated by the doctor before and after the placement of Breast Implants.

INTERFERENCE IN BREASTFEEDING – The presence of silicone in breast milk has been researched in several studies, with no significant levels of the substance found in mothers with silicone breast implants, when compared to women without implants. However, there is a possible interference in breastfeeding (reduced milk production) in women who have undergone surgery with breast implants. The frequency of interference with breastfeeding is higher in patients in whom the periareolar incision route was used.

INTERFERENCE WITH MAMMOGRAPHY; ULTRASOUND AND MAGNETIC RESONANCE IMAGING - Interference occurs when the presence of implants (or artifacts) causes an erroneous reading of the result.

SILIMED points out that the implant can interfere with the quality of the mammogram. Therefore, you must be guided and instructed to request the services of professionals who have experience with mammography techniques with implant, informing the radiologist of the need to adapt the mammographic compression and not forgetting to show the “SILIMED PRODUCT TRACEABILITY CARD - FOR PATIENT”.

Other methods, such as ultrasound and MRI, can be useful in conjunction with mammography because they do not require compression and allow examination from any angle.

According to recent studies, silicone breast implants reduce the image quality of echocardiographic visualizations. More data is needed, however SILIMED advises its patients to inform the technician responsible for the exam about the existence of silicone implants.

INTERFERENCE WITH SELF-EXAMINATION - Breast self-examination is a preventive technique used in an attempt to identify early stages of breast cancer. The doctor should guide you on how to distinguish the implant from the breast tissue during self-examination.

You should periodically self-examine for lumps, swelling, hardening, or changes in the shape of the implant, which may be signs of rupture. If any of these symptoms or persistent pain occurs, report it to your doctor.

It is important to note that the breast examination performed by the woman herself does not replace the physical examination performed by a health professional (doctor or nurse) qualified for this activity.

PRECAUTIONS - The risks of breast implant surgery may be greater if you have any of the following specific conditions. Tell your doctor if you have any of the following conditions:

- Autoimmune disease;
- Recent cancer history;
- Weakened immune system (users of immunosuppressive drugs);
- Conditions that interfere with blood clotting;
- History of severe allergy;
- Cardiovascular diseases;
- Diabetes;

Before surgery, you should have a detailed conversation with your doctor about your medical history.

CAUTIONS - In the first month after surgery, some activities can damage the implant and should be avoided, such as:

- Exposure to the sun;
- Sudden movements;
- Sports in general.

You should be aware that normal stress or trauma applied to the surgery site can cause the extrusion of the implant.

Wearing an appropriate bra and not doing violent exercises are minimum recommendations to be followed in the postoperative period.

Ask your doctor about activities he does not recommend in the postoperative process

You must inform the medical or pharmaceutical specialist that you have silicone implants before using topical medications (for example, steroids) in the breast area.

You need to continue to consult with your medical specialist to proceed with routine monitoring for breast cancer.

You should inform the doctor about the presence of the implant if any breast surgery is scheduled.

If you suspect any complications listed in this newsletter or any abnormal symptoms, contact your doctor.

DURABILITY - Like all breast implants, SILIMED Mammary Implants have a limited lifetime. This lifetime has not yet been determined by the scientific community. SILIMED establishes an average period of 10 years as the expected useful life, but this parameter can be changed when any reason that justifies it is raised.

The implant may need to be removed or replaced, which may require revision surgery.

Annual clinical follow-up of patients is recommended.

IMPORTANT FACTORS TO CONSIDER:

1. If you are going to undergo breast augmentation or reconstruction, be aware that breast implants are not considered to be lifelong devices and that breast implantation may not be a one-time surgery. It is possible that you may need one or more additional surgeries, in addition to medical monitoring throughout your life. It is also possible that you may need surgery to remove the implant, with or without replacement.

2. Many of the changes in your breast after implantation cannot be undone. If you choose to have your implant(s) removed in the future, you may experience unacceptable folds, drooping or wrinkling, loss of breast tissue or other alterations to your breast contour.

3. With breast implants, routine screening mammography will be more difficult and you will need to acquire additional images, which means more exposure to radiation.

4. Every SILIMED Breast Implant comes with a "SILIMED PRODUCT TRACEABILITY CARD - FOR PATIENT", which contains all information regarding the surgery. This card is for your safety, you must carry this card to facilitate medical care in an emergency.

This document is available on our website: www.silimed.com.br

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CLARIFICATION TERM FOR THE USE OF SILIMED MAMMARY IMPLANTS

Dear Patient,

In order to support the decision regarding the use of the SILIMED Mammary Implant, we have available on the company website the Information for Patients - SILIMED Mammary Implants, so that you can have access to important matters and clarify possible doubts with your doctor, before undergoing surgery.

Before and after surgery, you should consult with the doctor and undergo periodic examinations to monitor your health conditions.

It is important that you be aware of the need to correctly follow the guidelines provided by the doctor, to be successful in the surgery to place the SILIMED Mammary Implant.

Both the doctor and the anesthesiologist must have their consent to take the necessary measures to restore their clinical condition, if necessary.

The success of the surgery with the SILIMED Mammary Implant depends on the commitment to follow the medical recommendations, which will allow an increase in the projection or reconstruction of the breast and without major complications. Breast surgery can improve your self-esteem and quality of life.

Warnings:

- Breast implants have a limited life and may need to be removed or replaced, which may require revision surgery.

- Breast implants are associated with the development of a type of cancer of the immune system called anaplastic large cell lymphoma (BIA-ALCL), with greater frequency in textured breast implants; Breast implant patients are at increased risk of developing BIA-ALCL between the scar tissue and the fluid around the breast implant; BIA-ALCL usually takes several years to develop, but cases have been reported within one year after implantation; Typical symptoms to watch for include: swelling, breast tightness, pain, lumps or swelling of the breast, months or years after implantation surgery; The treatment of BIA-ALCL involves surgery to remove the implants and the surrounding scar tissue capsule. Based on the stage of cancer at diagnosis, some patients may require chemotherapy or radiation therapy. While some patients respond well to treatment, some die from BIA-ALCL.

- Some patients who received breast implants reported a variety of systemic symptoms including joint pain, fatigue, rash, memory loss, confusion, among others. Although the causes of these symptoms are unclear, some patients have reported relieving these symptoms by removing their implants and the surrounding scar tissue capsule

Confirmation

I, _____, identity _____, read the INFORMATION FOR PATIENT - SILIMED MAMMARY IMPLANTS made available by SILIMED, understood the product information and surgical procedures and had the opportunity to clarify with my physician any doubts about the nature of the product, its benefits, risks and potential complications, procedures and objectives of the surgery. I am aware of the possibility of the occurrence of adverse effects, such as those provided for in SILIMED Patient Information, and that the occurrence of adverse effects does not characterize a product defect, for warranty purposes.

Patient's Signature _____ date ____/____/____ ID Number: _____

Witness _____ date ____/____/____ ID Number: _____

Physician _____ date ____/____/____ Regional Board of Medicine: _____

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1st copy for the Patient / 2 nd copy for the Patient

CLARIFICATION TERM FOR THE USE OF SILIMED MAMMARY IMPLANTS

Dear Patient,

In order to support the decision regarding the use of the SILIMED Mammary Implant, we have available on the company website the Information for Patients - SILIMED Mammary Implants, so that you can have access to important matters and clarify possible doubts with your doctor, before undergoing surgery.

Before and after surgery, you should consult with the doctor and undergo periodic examinations to monitor your health conditions.

It is important that you be aware of the need to correctly follow the guidelines provided by the doctor, to be successful in the surgery to place the SILIMED Mammary Implant.

Both the doctor and the anesthesiologist must have their consent to take the necessary measures to restore their clinical condition, if necessary.

The success of the surgery with the SILIMED Mammary Implant depends on the commitment to follow the medical recommendations, which will allow an increase in the projection or reconstruction of the breast and without major complications. Breast surgery can improve your self-esteem and quality of life.

Warnings:

- Breast implants have a limited life and may need to be removed or replaced, which may require revision surgery.
- Breast implants are associated with the development of a type of cancer of the immune system called anaplastic large cell lymphoma (BIA-ALCL), with greater frequency in textured breast implants; Breast implant patients are at increased risk of developing BIA-ALCL between the scar tissue and the fluid around the breast implant; BIA-ALCL usually takes several years to develop, but cases have been reported within one year after implantation; Typical symptoms to watch for include: swelling, breast tightness, pain, lumps or swelling of the breast, months or years after implantation surgery; The treatment of BIA-ALCL involves surgery to remove the implants and the surrounding scar tissue capsule. Based on the stage of cancer at diagnosis, some patients may require chemotherapy or radiation therapy. While some patients respond well to treatment, some die from BIA-ALCL.
- Some patients who received breast implants reported a variety of systemic symptoms including joint pain, fatigue, rash, memory loss, confusion, among others. Although the causes of these symptoms are unclear, some patients have reported relieving these symptoms by removing their implants and the surrounding scar tissue capsule

Confirmation

I, _____, identity _____, read the INFORMATION FOR PATIENT - SILIMED MAMMARY IMPLANTS made available by SILIMED, understood the product information and surgical procedures and had the opportunity to clarify with my physician any doubts about the nature of the product, its benefits, risks and potential complications, procedures and objectives of the surgery. I am aware of the possibility of the occurrence of adverse effects, such as those provided for in SILIMED Patient Information, and that the occurrence of adverse effects does not characterize a product defect, for warranty purposes.

Patient's Signature _____ date ____/____/____ ID Number: _____

Witness _____ date ____/____/____ ID Number: _____

Physician _____ date ____/____/____ Regional Board of Medicine: _____

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1st copy for the Patient / **2nd copy for the Physician**

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