CE
1434

IUP 001-01-P-E R16

SILIMED 

IMPLANTE MAMÁRIO - GEL DE SILICONE - SUPERFÍCIE REVESTIDA COM ESPUMA DE POLIURETANO

PT

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INSTRUCTIONS FOR USE / INFORMATION ON THE PRODUCT

Review date: 30/12/2019

DESCRIPTION

MAMMARY IMPLANT – SILICONE GEL – POLYURETHANE FOAM COATED SURFACE

Consisting of a single, mechanically resistant silicone elastomer Shell coated with polyurethane foam. Contains a defined volume of transparent, high performance silicone gel (HSC)², developed so that its shape and consistency bring a natural breast look.

Low bleed treated Shells contribute to the reduction of silicone osmosis transudation through the elastomeric Shell barrier, where part of the Shell thickness is formed by an inner layer of elastomer especially impervious to silicone oils.

The polyurethane foam coated surface has scientific evidence¹ of its efficiency in greatly reducing the occurrence of capsular contracture. The polyurethane coating allows for better interaction between the implant and the surrounding tissue. The very uneven surface of the polyurethane foam disadvantages the formation of the linear fibrotic capsule in favor of the multiplanar formation of collagen fibrils. The formation of microcapsules around the polyurethane particles causes the contracture force to cease to have a single orientation, having multiple vectors, this force tends to cancel out, reducing capsular contracture.

It is presented in a set of types, which is the result of clinical experience gained from the medical profession.

Accessory: The Mammary Implant - Silicone Gel - Polyurethane Foam Coated Surface is supplied with an auxiliary glove for implantation. This accessory is not sold separately.

It is up to the physician to select the appropriate size and type of fixation to meet the clinical and aesthetic requirements of each case.

They are all supplied sterile and intended for single use only. The package should only be opened inside the operating room.

CONSTITUENT MATERIALS -The raw materials are medical grade and biocompatible. Environmental manufacturing conditions and production techniques are controlled by a strict quality management system.

- Shell: Elastomer, Composition of Polydimethylsiloxane and

¹ Pitanguy I, Salgado F, Radwanski HN & Stersa RM – Estágio Atual dos Implantes Mamários. Ver bras Cir. (Current Stage of the Mammary Implants. See bras. Cir.) 1991; 81(6): 291-299

Rebello C - Augmentation Mammoplasties With Polyurethane Foam-Coated Silicone Prosthesis. Rev da Soc Bras Cir Plast Est Reconstr, 1993; 6 (1,2,3): 47 - 57.
2 High Strength Silicone

Note 1: For more information regarding implants, see the table available on the cover of this Instruction for Use.

Note 2: Natural, Nuance and Enhance models may have a position indicator on the front and / or back to facilitate positioning of the product during the surgical procedure.

Note 3: The difference between the models is in the dimensions and / or volumes of the product.

Dimethyl Fluoro Silicone Copolymer, catalyzed by platinum compound

- Coating: Polyurethane foam adhered to the Shell by vulcanization process.

- Filler: Silicone gel (HSC)², catalyzed by a platinum complex.

Silicone may vary in color from opaque gray to translucent yellow.

Polyurethane foam may vary in color from cream to slightly brown.

PRESENTATION - Supplied sterile and pyrogen-free, in double packaging, packed in sealed outer box, where the documents related to the product are.

The outer packaging label contains, but is not limited to, the following information: product description, reference number, quantity of products per package, serial number, expiration date, product size and / or volume. The surface of the product is identified by the color used on its label as follows:

- Polyurethane Foam Coated Surface: Orange.

The symbols shown on the label are described on the form attached to this instruction for use.

STERILIZATION - One of the following sterilization methods is used:

- Ethylene oxide

The sterilization process and its expiration date are indicated on the packaging labels, meeting the requirements of the country to which the product is intended and which has granted registration for marketing. Each sterilization batch receives its individual confirmation.

The Mammary Implant is for single use only.

Forbidden to reprocess.

INTENDED USE / PURPOSE OF USE - SILIMED Silicone Gel Mammary Implants achieve their intended purpose in the following ways:

- in case of breast enlargement: increased volume and circumference of one or both breasts;

- in the case of breast reconstruction: return of the volume and circumference of one or both breasts removed after mastectomy or trauma or not developed due to any medical abnormality.

SILIMED Mammary Implant is a single-use product implantable through invasive surgery and its medical grade silicone elastomer composition and polyurethane foam lining, where applicable, come into contact with breast tissues for a prolonged use of the implant.

INDICATION OF USE - SILIMED Silicone Gel Mammary Implants are indicated for women for the following procedures:

- Breast enlargement. Breast augmentation includes preliminary breast augmentation to increase breast size, as well as revision surgery to correct or improve the outcome of preliminary breast augmentation surgery.

- Breast reconstruction. Breast reconstruction includes primary reconstruction to replace breast tissue that was removed due to cancer or trauma or that failed to develop properly due to

a serious breast abnormality. Breast reconstruction also includes revision surgery to correct or improve the results of a primary breast reconstruction surgery.

CONTRAINDICATIONS - The use of mammary implants or the implantation procedure is contraindicated in the following situations:

- Patients with infection anywhere in the body;
- Patients with cancer or premalignant condition not adequately treated;
- Pregnant or nursing patients;
- Patients with atopy;
- Patients with depressed immune system;
- Patients with local inflammation.

PRECAUTIONS - The physician should, in the preoperative evaluation of his patients, take into account: recent history of tumors in the implant region, in particular recurrent cancer or metastasis; autoimmune diseases; pathologies and treatments that affect the immune system or blood clotting and; history of severe allergy; diabetes and cardiovascular problems.

HANDLING AND USE INSTRUCTIONS

• Package Opening:

- 1 - Make sure that the external plastic cover has not been opened;
- 2 - Remove the double packaging from the inside of the sealed box. This should be carefully examined before use in the operating room. The product that has its packaging in any way violated should not be used;
- 3 - Separate the documents that accompany the product;
- 4 - Stick the adhesive labels with the product data to the hospital records, the doctor and also the file that will be delivered to the patient, as indicated by the symbols;
- 5 - Open the outer blister with access to the sterile inner blister containing the product. Be careful not to contaminate it outside the first one.
- 6 - Then open the sterile inner blister in the surgical field.

The high dielectric value of silicone can generate electrostatic charges, which are responsible for attracting particles in the environment, such as fluff and talc, to name a few. External contaminants that adhere to the Shell surface can generate reactions such as those caused by foreign bodies through increased fibrosis and fluid generation. For this reason, the precautions taken to open the package are very important.

SURGICAL PROCEDURE - Appropriate surgical techniques and procedures are necessarily the responsibility of the physician, who should evaluate them according to his own training and experience, observing updated and accepted techniques to minimize the occurrence of adverse reactions. The most common techniques and complications are described in this document. In addition, SILIMED has a scientific bibliography available from its representatives or directly from the head office.

It is also the doctor's responsibility to make sure that his staff is adequately trained to perform the required medical procedures.

The physician should discuss the potential risks and benefits of a preoperative mammogram and ultrasound for future reference.

In case of suspicious changes in these exams, we suggest investigation before or during implant placement.

Prior to the surgical procedure, the patient should be informed of all procedures to be submitted and all risks involved. The application of preoperative informed consent is recommended.

The physician should carefully evaluate the implant type, access route, pocket dissection and implant location, taking into account the patient's anatomy and the clinical and aesthetic requirements of each case.

Because they are filled with high cohesive silicone gel, SILIMED implants can be introduced through 3 access ways: inframammary, periareolar and transaxillary. The choice of route should be in accordance with the patient's anatomy and the surgeon's experience.

The inframammary route is technically the simplest due to its easy access to the subglandular, subfascial and submuscular planes. It is also easy to make the pocket, perform hemostasis, insert the implant and confirm its positioning.

The half-moon periareolar route requires a dissection bordering the mammary gland or sectioning it until it reaches the subglandular, subfascial and submuscular planes. Depending on the thickness of the gland and the surgeon's experience, it may be simpler or more laborious to make the pocket in the chosen plane, perform hemostasis and introduce and confirm implant placement. Creation of the asymmetrical pocket between the breasts can lead to a breast fold asymmetry. Periareolar introduction may also be associated with a higher incidence of infection and biofilm formation due to contact of the implant with the ruptured breast ducts during implant introduction, as well as increasing the risk of interference with breastfeeding.

Transaxillary access requires the dissection of a tunnel between the axillary incision and the mammary space. Therefore, the distance between the axillary access and the boundaries of the plane-independent pocket – subglandular, subfascial and submuscular - is greater. Some surgeons may find technical difficulty in performing hemostasis and accurately dissecting the boundaries of the mammary fold via this route. The introduction of the implant requires greater skill and care in its introduction in order to avoid axillary skin trauma and implant gel fracture due to excessive force. Incorrect implant placement can lead to implant asymmetry, rotation and displacement.

Regardless of the route chosen, the skin incision should be appropriate for correct visualization and dissection of the pocket. We suggest that the incision be between 5 and 7 cm, depending on the type of access chosen and the size of the implant. Disproportion between these aspects can lead to excessive force causing burns to the edges of the incision due to friction or fracture of the implant gel. Excessive pressure should be avoided during implantation.

Regarding positioning, three implant allocations can be used: retroglandular, retrofascial and retromuscular. All have advantages and disadvantages. These should be exposed and explained to the patient prior to implantation. As a general rule, the pre-pectoral (retroglandular and retrofascial) pathways preserve the pectoral musculature and have a lower incidence of postoperative pain. However, in patients undergoing the inclusion of large or small breast implants, the implants, especially their

edges, may be palpable and visible. The retromuscular plane has a thicker coverage, making it difficult to palpate and visualize the implant, but is more painful, changes muscle function and can lead to poor positioning of the implant with lateral and superior displacement.

Once the implants have been introduced, their placement should be confirmed. In round and conical implants, advance, maximum models, the symmetry of positioning between the breasts should be determined, avoiding postoperative misplacement, with lateral, cranial or caudal displacement of the implant, as well as groove asymmetries. In the case of natural, nuance and enhance models implants, besides confirming the symmetry of position between the breasts, as previously described, make use of the colored implant marking (visual indication) and embossed marking at the base of the implant (tactile indication), to determine the correct orientation, avoiding positions with the rotated or asymmetrical axis.

The plastic glove accessory aims to help the inclusion of the breast implant - surface coated with polyurethane foam in the implantation pocket. The implant must be allocated into the device through one of the openings and, at the physician's handling, the implant is introduced into the pocket. Once the implant has been placed in the desired position, this accessory should be removed.

The greatest immediate perioperative risks are known to be: hematoma, seroma and infection. Thus, correct asepsis, hemostasis and implant handling are recommended in order to reduce these complications. The use of drains is possible according to surgeon preference and / or patient need.

Basic postoperative recommendations include wearing a bra that is appropriate to the procedure performed, avoiding activities that may favor pain, bleeding or displacement of implants, and avoiding smoking.

For postoperative follow-up, patient and implant evaluation should be done through anamnesis and physical examination to detect possible changes in the implant. Imaging such as ultrasound and nuclear magnetic resonance are the best so far to assess the integrity of implants. We suggest performing an annual ultrasound for implant evaluation.

The physician should always have a spare product available during surgical procedures.

We suggest an annual clinical follow-up.

INSTRUCTIONS FOR REMOVAL OF THE IMPLANT - Breast implant removal is usually a simple procedure. However, possible previous complications may add some complexity to the surgery.

The access route for implant removal should be the most appropriate at the physician's discretion, based on previous examinations. After skin incision and dissection of the tissue planes, the plane used for implant placement should be accessed. Proceed with careful removal of the implant. Then close the layered planes and finalize the surgical procedure.

It is up to the physician to choose whether or not to remove the fibrotic capsule, assessing the impact of this withdrawal on patient morbidity in the postoperative period and also on possible interference with the interpretation of mammographic

examinations. In case of permanent withdrawal without further implant replacement, the physician should also evaluate the aesthetic result that the extraction of the capsule will cause.

In case of removal of the implant whose Shell suffered loss of integrity with contact of the gel with the patient's tissues, the physician should pay attention to the removal of the entire gel content and should also proceed with the removal of the fibrotic capsule.

Meticulous hemostasis during surgery is recommended in order to decrease the chances of postoperative hematoma and seroma. Excessive and persistent bleeding should be controlled throughout the procedure.

Damaged product must be returned to manufacturer for disposal.

IMPORTANT RECOMMENDATIONS / WARNINGS - Breast Implants may only be purchased by or on the order of a physician.

It is supplied sterile and, to maintain its aseptic condition during surgery, asepsis and cleanliness are required under the conditions of use.

The Implant and its packaging must be intact, otherwise the implant should not be used. It can only be used if found intact in its original form of manufacture, i.e. without any change in its original characteristics.

The implant may eventually have blisters or some whitening resulting from sterilization, which does not preclude its use.

Polyurethane foam is degraded, which consists of its mechanical and / or enzymatic fragmentation. The free slides remain contained in the fibrotic capsule. Until degradation is complete, the capsule is already properly formed without compromising implant integrity.

The implant surface should not be contaminated with talc, dust or oils. Before contacting the implant, the medical team should exchange the sterile surgical gloves in use for new pairs and then wash them with sterile, pyrogen-free saline.

Solutions containing iodine cannot contact the implant.

Marking solutions, such as methylene blue, used in the medical-surgical procedure must be properly sterilized to prevent biological contamination.

All perforating surgical instruments to be used in the region near the implant must be handled with care, so that it is not reached. Any accidental puncture, cut, or even scratching of the implant Shell will invalidate it for surgical use. It is not possible to repair damaged product. It must be replaced by another.

The physician should instruct his or her patient not to massage and / or vigorously manipulate the implant or to engage in intense physical activity for a specified period.

The patient should be advised to follow the most up-to-date medical recommendations regarding breast exams and to distinguish breast tissue implantation for best self-examination results.

Although definitive data are not yet available, the scientific literature route suggests that there is no direct cause-and-effect

correlation between breast implant placement for aesthetic augmentation and suicide frequency, and increased frequency is possibly correlated with intrinsic psychiatric morbidity. of these patients. Therefore, it is highly recommended that the patient's psychological stability be carefully evaluated by the physician before and after breast implant placement.

The Mammary Implant is for single use only. Forbidden to reprocess.

By law, it may not be reused or resterilized as it may damage the product and compromise its performance and safety.

Damaged product must be returned to manufacturer for disposal.

POTENTIAL COMPLICATIONS AND POSSIBLE SOLUTIONS

- Any patient undergoing surgery to introduce a foreign element into the body is subject to possible complications.

SILIMED delegates to physicians the responsibility to inform patients about the possibility of new surgery for implant removal or replacement, as well as possible adverse effects after implantation and appropriate treatment. At the same time, it invites the medical profession to report other findings in this regard.

SILIMED also delegates to the physician the adoption of the best solutions available for the complications of surgery. Among these, it is worth mentioning the meticulous hemostasis during the procedure, antibiotic therapy and removal or replacement of the implant.

PREDICTED COMPLICATIONS - For Mammary Implant, the potential complications of surgery, reported in the current medical-scientific literature, or directly to SILIMED, followed through the product risk management process - post-production phase, are as follows:

CHEST WALL FLATTENING - Few cases have been reported related to tissue expanders and occasionally to breast implants placed in the retromuscular position. It is associated with the surgical technique used. It may cause pain and injury to organs near the ribs and muscles of the thoracic region.

ALLERGENICITY (CUTANEOUS ERUPTION - RASH) - It may occur in the recent postoperative period and it is mainly characterized by tissue irritation, redness and rash. It is easily treatable and of transient duration.

CHANGE IN BREAST AND / OR NIPPLE SENSITIVITY - Uncommon. Corresponds to the increase (hyperesthesia) or decrease (hypoesthesia) of sensitivity of the implanted region due to the surgical incision. This is usually a reversible process over time, but it can last two years or more.

TISSUE ATROPHY - Associated with large volume breast implants, when inserted in pockets disproportionate to the volume of these implants. Mechanical forces put pressure on the parenchyma and subcutaneous tissue and may cause atrophy of these tissues over time.

CALCIFICATION - It is characterized by the presence of insoluble calcium (and / or magnesium) salt deposits, especially calcium carbonate and phosphate in the fibrotic capsule. May cause pain and hardening of healing tissue. Of unknown origin, it has been observed in tissues surrounding the implant,

sometimes requiring implant removal.

CYSTS / NODULES - Breast nodules and cysts are changes that occur in the breast tissue. Breast lumps are breast protuberances. They can usually be caused by fibrocystic changes, fibroadenomas, cysts and cancer. They may also occur due to implant rupture if extravasation occurs and the silicone gel migrates to lymph nodes or other parts of the body. The cyst is characterized as a closed sac / cavity filled with fluid lined with an epithelium.

IATROGENIC COMPLICATIONS - Certain situations such as folds and perforations in the breast implant, among others, may be caused by the use of improper surgical technique.

CAPSULAR / FIBROSIS CONTRACTURE - It occurs when the healing tissue or fibrotic capsule becomes thicker and more retractable and compresses the implant. The fibrotic capsule forms around the implant due to the body's reaction to the presence of a foreign body. Factors such as infection, hematoma and seroma may favor the appearance of capsular contracture. The probability of occurrence decreases with meticulous pocket hemostasis and care in implant handling. SILIMED advises against the external (closed) capsulotomy procedure, which may cause implant rupture.

SUTURE DEHISCENCE - Suture dehiscence is the separation of the layers of the surgical wound. It may be partial and superficial or complete when there is complete disruption of all layers of the tissue. Usually this effect begins with a small seroma. In most cases, this is a problem related to the surgical technique and / or use of a larger than appropriate implant.

IMPLANT DISPLACEMENT / BAD POSITIONING - It may occur with patient discomfort and distortion as desired, which is a problem linked to the operative technique, ineffective formation of the fibrotic capsule or trauma. The physician should evaluate the need for implant removal and / or likely replacement.

PAIN - Some patients may experience pain in the postoperative period, which should be investigated immediately. It may occur mild to severe and short to long term. It may result from infection, inflammation, seroma, calcification and / or improper surgical technique.

WRINKLE - Breast implant surface ripple that can be seen and / or felt.

EROSION - After implantation, tissue erosion may occur near the implant, caused by device folds and other deformities that lead to over compression of the tissue.

EXTRUSION - Unstable or compromised cover tissue, excessive stress, or trauma occurring at the implant site and / or disruption of the healing process may result in implant extrusion. It has been reported as an uncommon complication in both patients undergoing augmentation and reconstruction surgery.

GALACTORREIA - May be the result of hormonal disorders or the action of medications. It causes discomfort and embarrassment, affecting the emotional state of the patient.

GRANULOMA - It is associated with the formation of granular tissue due to the presence of various cell types such as macrophages, lymphocytes and fibroblasts. It may lead to

a chronic inflammatory lesion, tumor and consequent surgical intervention for implant removal.

HEMATOMA - Extravasation of blood and its confinement in an organ or tissue (hematoma) can occur in any type of surgery. In its occurrence, it should be carefully evaluated to define the conduct to be adopted. At risk of being a precursor of infection and increased fibrosis, if not resorbed, the bruise should be removed. The main preventive measure is perfect hemostasis during surgery and postoperative rest.

INFECTION (BIOLOGICAL CONTAMINATION) - Uncommon and unlikely due to implantation if the recommendations for use and asepsis associated with implantation surgery are observed. Includes the possibility of secondary necrosis and implant extrusion. Local culture is recommended, followed by antibiotic use, drainage and eventually implant removal.

INTERFERENCE WITH MAMMOGRAPHY; ULTRASOUND AND MAGNETIC RESONANCE - Interference occurs when the presence of implants (or artifacts) causes an erroneous reading of the result, for example: false positive diagnosis of rupture, resulting in unnecessary removal. Breast implants can interfere with the diagnosis of breast cancer by mammography, covering up any suspicious stains in the breast tissue, masking the results and delaying the initiation of anti-cancer treatment.

BREASTFEEDING INTERFERENCE - The presence of silicone in breast milk has been researched in several studies and no significant levels of breastfeeding have been found in mothers with silicone breast implants when compared to women without implants. However, there is a possible interference with breastfeeding (reduction in milk production) in women who underwent breast implants. The frequency of interference with breastfeeding is higher in patients whose periareolar breast incision route was used.

NECROSIS - Resulting from excessive skin tension caused by the presence of the implant or surgical trauma. The physician should evaluate the need for implant removal and / or likely replacement.

LOSS OF MECHANICAL INTEGRITY / PRODUCT FAILURE - In the event of loss of mechanical integrity of the breast implant, the patient may undergo surgery for replacement depending on the degree of danger. Among these situations, we can mention the implant detachment, the spontaneous rupture of its membrane, among others.

PNEUMOTHORAX - This is the collapse of the lung due to a sudden change in pressure within the chest cavity. In the case of breast implants, this may occur when the implantation shop is formed due to surgical trauma, especially in the retromuscular space. It has as consequences: tachycardia, chest pain, cough, shortness of air, cyanosis, agitation, among others.

INFLAMMATORY RESPONSE - The presence of an element foreign to the body triggers an inflammatory process, which may be more or less intense. Intense inflammation may also be due to inappropriate trauma or procedures. Causes pain, redness, local heat and swelling. This effect can be treated pharmacologically.

UNSATISFACTORY AESTHETIC RESULT - It is associated with patient and / or physician dissatisfaction with the outcome of

breast implant placement surgery. Some unsatisfactory possible aesthetic results are: ptosis, sagging, asymmetry, folds, folds, implant wrinkles, stretch marks, adhesion and hypertrophic scars. The physician should evaluate the best action to be taken, and there is no need for surgical intervention to replace the implant.

Rupture - It can be caused by the natural wear of the implant during implantation time, damage by a sharp instrument, excessive exertion during surgery, trauma or extreme physical pressure, compression during mammography or severe capsular contracture. Most of them occur during surgery. In case of suspected implant rupture in an asymptomatic patient ("silent rupture"), a careful evaluation should be made to confirm the situation and decide on the procedures to be taken, including the possibility of implant replacement.

SEROMA - Accompanied by edema and pain, which may subside with rest, immobilization and application of compresses, eventually drainage and even removal of the implant. The seroma is not always acute. There may be a minimal, continuous, undetected seroma that appears months later as "late seroma". There are even "late seromas" of unknown origin.

SYMPTOMS ASSOCIATED WITH MACROMASTY - They are associated with large breasts. Some of the symptoms are: backache, pain in breasts, hands (including numbness), neck and headache; skin irritation and shoulder markings. Breast implant volume should be considered with respect to the amount of breast tissue present, breast shape and the patient's desire. The methodologies used to estimate whether the volume of breast implants is adequate to the needs and expectations of patients are of great importance in choosing the breast implant.

GEL BLEEDING - This is associated with a silicone gel-filled breast implant. It is the process of diffusion and migration of the gel through the whole membrane (transudation) or rupture of the implant to adjacent tissues. It may cause local complications such as: reduction in size, asymmetry, change in breast shape, pain, capsular contracture or formation of nodules and granulomas. SILIMED high performance silicone gel minimizes the occurrence of this adverse effect.

THROMBOSIS - Uncommon, usually associated with the surgical procedure or patients with obesity, malignancies, varicose veins, traumatic or iatrogenic venous lesions, patients bedridden or immobilized by a systemic or local disease in the limb, patients with heart disease, extrinsic vein compressions, pregnancy and puerperium, use of hormonal contraceptives and hematologic diseases. May cause local complications, and implant removal is required.

LARGE CELL ANAPLASTIC LYMPHOMA ASSOCIATED WITH BREAST IMPLANT (BIA-ALCL) - BIA-ALCL is a type of non-Hodgkin's lymphoma, which is cancer of the immune system (not breast cancer), most often found in the tissue. Scarring that surrounds the implant, also called a capsule, or in the fluid that is formed between the implant and the capsule (seroma). The occurrence of late seromas is widely considered by the plastic surgery community to be associated with the pathogenesis of ALCL.

The pathogenesis of BIA-ALCL is still unknown and its

development is slow. The mechanism believed to be related to the development of BIA-ALCL is the exaggerated inflammatory response of the body (foreign body reaction). Another method that can trigger BIA-ALCL is bacterial biofilm infection around breast implants, leading to inflammatory response. An important factor that may be associated with the development of the disease is the genetic issue.

Although an increase in the prevalence of BIA-ALCL has been observed, its incidence is currently low. However, although any type of breast implant, regardless of surface, filling or manufacturer, may be associated with the development of ALCL, the prevalence is higher in textured surface implants.

Different treatments have been reported, such as implant removal and capsulotomy, lymph node dissection, chemotherapy, radiotherapy and even autologous stem cell transplantation. However, there is no universal treatment. If diagnosed, individual treatment should be performed. Women should monitor their breast implants and contact their doctors if they notice any changes, and those considering breast implant surgery should discuss the risks and benefits with their doctors.

Other Complications: Emotional Alteration, Discomfort, Numbness, Edema, Erythema, Breast Hardening, Bleeding, Swollen Lymph Nodes (Axillary Lymph Nodes) also Called Lymphadenopathy, Tissue Infiltration, Irritation, Tissue Injury, Palpability, Itching, Flushing, Skin Tension, Implant Visibility, Asymmetry, Ptosis, Late Healing, Upper Pole Filling.

Some complications have been reported in the medical-scientific literature, but so far there is no scientific evidence of a possible cause-and-effect relationship between these events and silicone implants. Some examples are:

BREAST CANCER - According to current medical-scientific literature, women undergoing silicone implant surgery are not at higher risk of developing breast cancer.

IMMUNOLOGICAL RESPONSE - In studies conducted so far, no evidence has been found that placement of silicone implants leads to connective or autoimmune diseases. Current literature suggests a possible link between implant rupture and fibromyalgia and connective tissue disease. Further studies are needed.

TERATOGENICITY OR OTHER REPRODUCTIVE EFFECTS- It is the ability to cause birth defects in the fetus. The review of the medical and scientific literature indicates that studies have not shown any evidence of teratogenicity or other reproductive effects associated with breast implants.

EFFECTS ON CHILDREN WHO ARE BORN FROM MOTHERS WITH MAMMARY IMPLANTS - Recent studies have shown that the amount of silicone in milk found in women who have silicone implants is not significantly different from women who do not have silicone implants. However, the review of the medical and scientific literature shows no evidence that silicone gel breast implants can have any harmful effects on the children of implanted women.

SUICIDE - Epidemiological studies indicate that women who have undergone cosmetic augmentation mammoplasty have a higher suicide rate than the general population. Although

there are no definitive data yet, the scientific literature route suggests that there is no direct cause and effect correlation between breast implant placement for aesthetic augmentation and suicide frequency, and the increase in this frequency is possibly correlated with intrinsic psychiatric morbidity of the patients. Therefore, it is highly recommended that the patient's psychological stability be carefully evaluated by the physician before and after breast implant placement.

DURABILITY - Like all breast implants, SILIMED Breast Implants have a limited lifespan. This lifetime has not yet been determined by the scientific community. SILIMED establishes a Medium period of 10 years as the expected useful life, but this parameter can be changed when any justifying reason is raised.

The implant may need to be removed or replaced, which may require revision surgery.

Annual clinical follow-up of patients is recommended.

SILIMED PATIENT INFORMATION -The following information is intended for all patients who intend to undergo augmentation or reconstructive surgery using silicone breast implants. It provides information on key questions surrounding the issue, including potential complications, benefits and risks, activities that may damage the implant, the possible need for device removal / replacement, and the expected durability of the breast implant.

The physician responsible for the procedure should read this section carefully and is responsible for resolving all patient questions prior to the SILIMED mammary implant procedure.

Annual clinical follow-up of these patients is recommended. If you suspect any complication listed in this newsletter or any abnormal symptoms, contact your doctor.

EXPECTED BENEFITS - A breast implant surgery can bring great benefit by providing successful breast reconstruction, augmentation and surgical revision. Breast reconstruction / augmentation also results in improvements in quality of life. In addition, studies show that breast implant reconstruction has been a help in breast cancer recovery, as well as reducing emotional stress by helping the body to look more natural as opposed to not having surgery reconstructive surgery or use external prosthesis.

PATIENT CLARIFICATION AND CONSENT - Considering the risks inherent in surgery, with or without the use of implants, and the possible complications related to them, SILIMED relies on physicians to clarify their patients regarding the existing risk-benefit balance.

Each patient should receive SILIMED Patient Information during their initial consultation to allow sufficient time prior to surgery to adequately read and understand the important risk information, follow-up recommendations, and benefits associated with breast implant placement surgeries. To successfully document patient clarification, the SILIMED Breast Implant Clarification Term must be signed by the patient and physician and attached to the patient's record.

INTERFERENCE WITH MAMMOGRAPHY; UNTRASOUND AND MAGNETIC RESSONANCE - Interference occurs when the presence of implants (or artifacts) causes an erroneous reading

of the result, for example: false positive diagnosis of rupture, resulting in unnecessary withdrawals.

SILIMED points out that the implant may interfere with the quality of mammography. Therefore, the patient should be oriented and instructed to request the services of professionals who have experience with implanted mammography techniques, informing the radiologist about the need to adapt the mammography compression and not forgetting to show SILIMED "PRODUCT TRACEABILITY CARD" - FOR PATIENT".

Other methods such as ultrasound and magnetic resonance imaging may be useful in conjunction with mammography because they do not require compression and allow examination of any angle.

According to recent studies, silicone breast implants reduce the image quality of echocardiographic views. More data is needed, however SILIMED advises its patients to inform the technician responsible for the existence of silicone implants.

PRECAUTIONS - In the first month after surgery, some activities may damage the implant and should be avoided, such as:

- Sun exposure;
- Sudden movements;
- Sports in general.

You should be aware that normal stress or trauma applied at the surgical site may cause implant extrusion.

The use of proper bra and not doing violent exercises are minimum recommendations to be followed in the postoperative period.

The patient should inform the medical or pharmacist specialist who has silicone implants before using topical breast medicine.

The patient needs to continue to consult with the medical specialist to proceed with routine monitoring to detect breast cancer.

The patient should inform the doctor of the presence of the implant if any breast surgery is scheduled.

If you suspect any complication listed in this newsletter or any abnormal symptoms, especially if trauma or compression occurs (e.g. excessive breast massage through some sports or seat belt use), contact your doctor.

IMPORTANT FACTORS TO CONSIDER -

1. If you are undergoing breast augmentation or reconstruction, be aware that breast implants are not considered lifelong devices and that breast implantation may not be a single surgery. You may need one or more additional surgeries and have to schedule doctor visits throughout your life. You may also need surgery to remove the implant, with or without replacement.
2. Many of the changes to your breast after implantation cannot be undone. If you subsequently choose to have the implant (s) removed, you may experience unacceptable plecting, unacceptable wrinkling, unacceptable wrinkling, loss of breast tissue, or other changes in breast contour.
3. With breast implants, routine evaluation mammography will be more difficult and you will need additional images,

which means more radiation and more time.

4. Every SILIMED Breast Implant comes with a "Silimed Patient Traceability Card - For Patient" containing all information regarding surgery. This card is for your safety, you should carry this card to facilitate medical attention in an emergency.

STORAGE / CONSERVATION / TRANSPORTATION - SILIMED products and packaging are durable and, if pocketed at room temperature and according to the directions printed on the packaging, will not be damaged. Because they are sterile medical products, it is recommended that they are not transported or pocketed with other types of materials that could cause physical or chemical damage to the product, which would prevent their use.

WARRANTY - SILIMED replaces any product in its line that has a proven manufacturing defect, provided it is properly identified by the purchaser and after examination in its laboratory.

The adverse effects presented here, the inappropriate surgical use and the color variation mentioned in the item "Constituent Materials", are not considered product defects.

The SILIMED product warranty does not cover the simple decision of the patient or physician to change the product.

TRACEABILITY - SILIMED maintains accurate, unit-by-unit records of raw materials, manufacturing steps, environmental conditions, operating conditions, and quality testing by assigning a unique number to each product that fully identifies it at any time. In this regard, it is absolutely necessary that any claim be forwarded with the SERIAL NUMBER (SN) or BATCH NUMBER, where applicable, of the unit claimed on the packaging.

SILIMED provides adhesive labels with product data that must be affixed to the hospital and physician records as indicated by the symbols.

Also provided are Patient Registration Form, to be signed by the patient and returned to SILIMED by the physician, and SILIMED Product Traceability Card, to be delivered by the physician to the patient. These documents contain product data, which allows traceability.

RETURN OF PRODUCTS - Products contaminated with biological material must be properly decontaminated by the physician or health care facility prior to shipment and returned to SILIMED following applicable legislation.

No products must be returned without prior authorization from SILIMED.

Dear Doctor: In order to facilitate understanding when handling our products, in this printed notice you will find international symbols, together with those developed at SILIMED, that appear on our labels, packaging, instructions for use and/or catalogues, as well as their meanings.

	Single use		Do not use if package is damaged
	Lot number		Volume of the product
	Date of manufacture		Consult instructions for use
	Ethyleneoxide sterilization		Keep away from sunlight
	Steam or dry-heat sterilization		Dimensions of product
	Manufacturer		CE mark
	Use by		Quantity of product
	Serial number		ANVISA registration number
	Sterile		Label for tracing to the doctor's dossier
	Reference number		Verified by the QC

	Caution, consult accompanying documents		INMETRO seal for identification of conformity
	Authorised representative in the european community		Label for tracing to the hospital dossier
	Label for tracing the patient's dossier		Medical device
	Keep dry		Temperature limit
	This way up		Humidity limit
	Right Side		Unique Device Identifier
	Left Side		Portuguese
	Spanish		Italian
	Russian		English

EN

INFORMATION FOR THE PATIENT - SILIMED MAMMARY IMPLANTS

Silimed has developed this leaflet to help you understand the benefits and risks related to augmentation or reconstructive surgery using silicone breast implants. Once you have read this leaflet carefully, talk to your doctor before making any decisions. If you decide on the procedure, you and your doctor should sign the Clarification Form prior to surgery, confirming that you have read and understood the information provided about the benefits and risks of silicone breast implant surgery. The SILIMED Breast Implant Clarification is at the end of this leaflet.

DESCRIPTION - SILIMED Mammary Implants consist of a single outer elastomer membrane filled with medical grade silicone gel; they may have a smooth, textured surface or even coated with polyurethane foam. They are all supplied sterile and intended for single use only.

These implants are available in various sizes and shapes to help each woman achieve the best result for her body. The doctor should guide you in choosing the correct implant to achieve the desired aesthetic effect.

PREDICTED BENEFITS – Breast implant augmentation / reconstruction surgery can be of great benefit by providing high-satisfaction breast reconstruction, augmentation, and surgical revision. Patients choose primary breast augmentation surgery to increase the size and proportion of their breast (s). In addition, patients choose revision surgery (replacement of an existing breast implant) to correct or improve the results of a primary augmentation surgery. According to the medical and scientific literature, most women who underwent breast implant surgery reported high levels of satisfaction with their body image, shape and size. Additionally, after surgery, patients reported higher self-esteem, and felt more feminine and attractive. Breast implant reconstruction surgery also results in improved quality of life. In addition, studies show that breast reconstruction with breast implantation has been an aid in breast cancer recovery.

PATIENT CLARIFICATION AND CONSENT - Considering the risks inherent in surgery, with or without the use of implants, and the possible complications related to these, SILIMED relies on physicians to clarify their patients regarding the existing risk-benefit balance. Each patient should receive SILIMED Mammary Implant Patient Information during their initial consultation to allow sufficient time prior to surgery to properly read and understand the important risk information, follow-up recommendations, and benefits associated with breast implant placement surgeries. To successfully document patient clarification, the SILIMED Mammary Implant Use Clarification Form must be signed by the patient and physician and attached to the patient's record.

PREDICTED COMPLICATIONS - Breast implant surgery offers great benefits. However, any type of surgery involves risks. There are some local complications (problems near the surgery site) that may occur after surgery. Silimed delegates to physicians the responsibility of informing patients about possible adverse effects, consequences and appropriate treatment. The following are the potential risks and complications reported in the literature associated with breast implant surgery of various models and manufacturers, as well as the possible effects related to these risks.

Breast and / or nipple sensitivity alteration, Chest-wall flattening, allergenicity (rash), emotional alteration / depression, tissue atrophy, calcification, cysts / nodules / tumor, iatrogenic complications, capsular contracture / fibrosis, dehiscence suture, discomfort, implant dislocation, pain, numbness, edema (swelling), erythema, (hyperemia), hardening, wrinkling, erosion, extrusion, galactorrhea, granuloma, hematoma, hemorrhage, implant removal, infection (biological contamination) and its symptoms, swollen lymph nodes (axillary lymph nodes), tissue infiltration, irritation, tissue injury, breastfeeding interference, poor positioning, silicon migration (lymphadenopathy and silicoma), necrosis, loss of mechanical integrity / product failure, palpability, pneumothorax, pruritus, local inflammatory response, revision surgery, poor esthetic outcome (such as ptosis, sagging, asymmetry, folds, folds, adhesion, stretch marks and hypertrophic scars), flushing, rupture (including cases of "silent" rupture), secretion, seroma / fluid, symmastia, symptoms associated with hypertrophy, skin tension, transudation of the gel bleeding, thrombosis, vesicles, visibility of the implant, filling of the upper pole.

LARGE CELL ANAPLASTIC LYMPHOMA ASSOCIATED WITH BREAST IMPLANT (BIA-ALCL) - BIA-ALCL is a type of non-Hodgkin's lymphoma, which is cancer of the immune system (not breast cancer), most often found in the tissue scarring that surrounds the implant, also called a capsule, or in the fluid that forms between the implant and the capsule (seroma). The occurrence of late seromas is widely considered by the plastic surgery community to be associated with the pathogenesis of ALCL.

The pathogenesis of BIA-ALCL is still unknown and its development is slow. The mechanism believed to be related to the development of BIA-ALCL is the exaggerated inflammatory response of the body (foreign body reaction). Another pathway that can trigger BIA-ALCL is bacterial biofilm infection around breast implants, leading to inflammatory response. An important factor that may be associated with the development of the disease is the genetic issue.

Although an increase in the prevalence of BIA-ALCL has been observed, its incidence is currently Low. However, although any type of breast implant, regardless of surface, filling or manufacturer, may be associated with the development of ALCL, the prevalence is

higher in textured surface implants.

Different treatments have been reported, such as implant removal and capsulotomy, lymph node dissection, chemotherapy, radiotherapy and even autologous stem cell transplantation. However, there is no universal treatment. If diagnosed, individual treatment should be performed. Women should monitor their breast implants and contact their doctors if they notice any changes, and those considering breast implant surgery should discuss the risks and benefits with their doctors.

Other complications have been reported in the medical-scientific literature, but so far, there is no scientific evidence of a possible cause and effect relationship between these events and silicone implants. Some examples are:

BREAST CANCER - According to current medical-scientific literature, women undergoing silicone implant surgery are not at higher risk of developing breast cancer.

TERATOGENICITY OR OTHER REPRODUCTIVE EFFECTS - It is the ability to cause birth defects in the fetus. The review of the medical and scientific literature indicates that studies have not shown any evidence of teratogenicity or other reproductive effects associated with breast implants.

EFFECTS ON CHILDREN BORN OF MOTHERS WITH MAMMARY IMPLANTS - Recent studies have shown that the amount of silicone in milk found in women who have silicone implants is not significantly different from women who do not have silicone implants. However, the review of the medical literature does not show evidence that silicone gel breast implants can have any harmful effects on the children of implanted women.

SUICIDE - Epidemiological studies indicate that women who have undergone cosmetic augmentation mammoplasty have a higher suicide rate than the general population. Although there are no definitive data yet, the scientific literature route suggests that there is no direct cause and effect correlation between breast implant placement for aesthetic augmentation and suicide frequency, and the increase in this frequency is possibly correlated with intrinsic psychiatric morbidity of these patients. Therefore, it is highly recommended that the patient's psychological stability be carefully evaluated by the physician before and after breast implant placement.

BREASTFEEDING INTERFERENCE - The presence of silicone in breast milk has been researched in several studies and no significant levels of breastfeeding have been found in mothers with silicone breast implants when compared to women without implants. However, there is a possible interference with breastfeeding (reduction in milk production) in women who underwent surgery with breast implants. The frequency of interference with breastfeeding is higher in patients whose periareolar breast incision route was used.

INTERFERENCE WITH MAMMOGRAPHY; ULTRASOUND AND MAGNETIC RESONANCE - Interference occurs when the presence of implants (or artifacts) causes an erroneous reading of the result, for example: false positive diagnosis of rupture, resulting in unnecessary withdrawals.

SILIMED points out that the implant may interfere with the quality of mammography. Therefore, the patient should be oriented and instructed to request the services of professionals who have experience with implanted mammography techniques, informing the radiologist about the need to adapt the mammography compression and not forgetting to show the "PRODUCT TRACEABILITY CARD" SILIMED - FOR PATIENT".

Other methods such as ultrasound and magnetic resonance imaging may be useful in conjunction with mammography because they do not require compression and allow examination of any angle.

According to recent studies, silicone breast implants reduce the image quality of echocardiographic views. More data is needed, however SILIMED advises its patients to inform the technician responsible of the existence of silicone implants.

INTERFERENCE WITH SELF-EXAMINATION - Breast self-examination is a prevention technique used in an attempt to identify early stages of breast cancer.

The doctor should advise you on how to distinguish breast tissue implant during self-examination. You should periodically self-examine for lump tracking, swelling, hardening, or change in implant shape, which may be signs of rupture. If any of these symptoms or persistent pain occurs, report it to your doctor.

It is important to note that the woman's own breast exam does not replace the physical exam performed by a health professional (doctor or nurse) qualified for this activity.

PRECAUTIONS - The risks of breast implant surgery may be greater if you have any of the following specific conditions. Tell your doctor if you have any of the following conditions:

- Autoimmune disease;
- Recent history of cancer;

- Weakened immune system (users of immunosuppressive drugs);
- Conditions that interfere with blood clotting;
- History of severe allergy;
- Cardiovascular diseases;
- Diabetes;

Before surgery, you should have a detailed conversation with your doctor about your medical history.

CAUTIONS -In the first month after surgery, some activities may damage the implant and should be avoided, such as:

- Sun exposure;
- Sudden movements;
- Sports in general.

You should be aware that normal stress or trauma applied at the surgery site may cause implant extrusion.

The use of proper bra and not doing violent exercises are minimum recommendations to be followed in the postoperative period.

Ask your doctor about activities he does not recommend postoperatively.

The patient must inform the medical specialist or pharmacist if they are bearers of silicone implants before they use topical medications on the breast region.

You need to continue to consult with the medical specialist to proceed with routine monitoring to detect breast cancer.

You should inform the doctor about the presence of the implant if any breast surgery is scheduled.

If you suspect any complications listed in this newsletter or any abnormal symptoms, contact your doctor.

DURABILITY - Like all breast implants, SILIMED Mammary Implants have a limited shelf life. This lifetime has not yet been determined by the scientific community. SILIMED establishes a Medium period of 10 years as the expected useful life, but this parameter can be changed when any justifying reason is raised.

The implant may need to be removed or replaced, which may require revision surgery.

Annual clinical follow-up of patients is recommended.

IMPORTANT FACTORS TO CONSIDER:

1. If you are undergoing breast augmentation or reconstruction, be aware that breast implants are not considered lifelong devices and that breast implantation may not be a single surgery. You may need one or more additional surgeries and medical follow-up throughout your life. You may also need surgery to remove the implant, with or without replacement.

2. Many of the changes to your breast after implantation cannot be undone. If you choose to have your implant(s) removed in the future, you may experience unacceptable folds, drooping or wrinkling, loss of breast tissue or other alterations to your breast contour.

3. With breast implants, routine screening mammography will be more difficult and you will need additional imaging, which means more radiation exposure.

4. Every SILIMED Mammary Implant comes with a "SILIMED PATIENT PRODUCT TRACEABILITY CARD" which contains all information regarding the surgery. This card is for your safety, you should carry this card to facilitate medical attention in an emergency.

This document is available on our website: www.silimed.com.br

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1st copy for the Patient / 2nd copy for the Physician

CLARIFICATION TERM FOR THE USE OF SILIMED MAMMARY IMPLANTS

Dear Patient,

To support the decision regarding the use of the SILIMED Mammary Implant, we have available on the company website the Patient Information - SILIMED Mammary Implants, so that you can have access to important matters and clarify possible doubts with your doctor, before undergoing surgery.

Before and after surgery, you should consult your doctor and have periodic examinations to monitor your health.

It is important that you be aware of the need to properly follow the guidelines provided by your physician to be successful in SILIMED Breast Implant placement surgery. Breast implants have a limited shelf life and may need to be removed or replaced, which may require revision surgery.

Both the physician and anesthesiologist must have their consent to take the necessary steps to re-establish their clinical condition if necessary.

The success of surgery with SILIMED Mammary Implant depends on the commitment to follow the medical recommendations, which will allow the breast to be enlarged or reconstructed without major complications. Breast surgery can improve your self-esteem and quality of life.

Confirmation

I, _____, identity at _____, read the INFORMATION FOR PATIENT - SILIMED MAMMARY IMPLANTS made available by SILIMED, understood the product information and surgery procedures and had the opportunity to clarify with my doctor any doubts about the nature of the product, its benefits, risks and potential complications, procedures and objectives of the surgery. I am aware of the possibility of adverse effects, such as those provided in the SILIMED Patient Information, and that the occurrence of adverse effects is not a product defect for warranty purposes.

Patient's Signature _____ date ____/____/____ ID Number: _____

Witness _____ date ____/____/____ ID Number: _____

Physician _____ date ____/____/____

Regional Board of Medicine: _____ Signature: _____

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1st copy for the Patient / **2nd copy for the Physician**

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